

D-Lux+

Dual-Wavelength Cordless Curing Light

**BUY 1, GET 1 FREE OR
GET 20 DIAFIL/40 DIAFIL FLOW FREE***

Scan for video



- 1 Powerful 2,400 mW/cm² output
- 2 Dual wavelength: 385 - 515 nm
- 3 360° swivel light probe
- 4 5 versatile curing modes
- 5 Built-in Light Intensity Sensor
- 6 High capacity 2,600 mAh battery
- 7 Wireless inductive charger

[4008-1001] D-Lux+ Regular Kit

Contains: D-Lux+ hp & charger, light probe and protector, disposable sheaths (200/box), power supply (adapter)

*DiaFil syringe refills or DiaFil Flow syringe refills. Choices of shades.



DiaFil™

Light-Cured Nanohybrid Composite

**BUY 2, GET 1 FREE OR
BUY 4, GET 4 FREE!**



*of the same value

- 1 Nanohybrid composite
- 2 Easy to place/sculpt and does not slump
- 3 Does not stick to instruments
- 4 Adapts well to margins
- 5 Wear & high fracture resistant
- 6 12 shades with remarkable color stability
- 7 High tensile strength & low shrinkage
- 8 Durable restoration w/ long-lasting gloss shine

SYRINGE REFILL (1 x 4g)

Shades: A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, A20, A30



CAPSULE REFILL (20 x 0.25g)

Shades: A1, A2, A3, A3.5, B2

Dia-X™ Bond UNIVERSAL

Single Component Universal Bonding Agent

BUY 1, GET 1 FREE!

Scan for peer review



- 1 Versatile with all etching techniques
- 2 All direct and indirect restorations
- 3 Light-, self- & dual-cure compatible
- 4 Virtually no post-op sensitivity
- 5 Advanced 10-MDP monomer
- 6 Convenient: true 1-bottle bond
- 7 Uncompromising bond strength
- 8 Stable at room temperature

[A2001-2201] Dia-X Bond Universal

Contains: 1 x 5 ml bottle



DiaFil Flow™

Light-Cured Flowable Composite

**BUY 2, GET 1 FREE OR
BUY 4, GET 4 FREE!**

- Excellent color adaptability
- Pit & fissure sealant
- Excellent tensile strength and low shrinkage
- Highly radiopaque

REFILL PACK (2g syringe & 10 tips)

Shades: A1, A2, A3, A3.5, B2, B3, C2, A20, A30.



DiaPlus™

One-Step Bonding Agent (5th gen)

**Buy 1 bottle of DiaPlus,
Get 1 DiaEtch Reg. Pack FREE!**

- Designed for both direct and indirect applications (enamel, dentin, porcelain, ceramics, amalgam, etc.)
- Low film thickness prevents adhesive show-through
- Ethanol-based bonding agent: exceptional bond strength
- Compatible w/ dual-/self-cure composites & core materials

[A2001-2101] DiaPlus Bonding Agent

Contains: 1 x 5ml bottle



DUO-GUN

Cordless Backfill Obturation Device



Buy & Get:
1 bottle GP Obturators
1 DiaPrep Pro (EDTA)



DIA-DUO

Complete Cordless Obturation System



Buy & Get:
1 bottle GP Obturators
1 box of Duo-Gun tips
1 DiaPrep Pro (EDTA)



DUO-PEN

Cordless Warm Vertical Compaction Device



Buy & Get:
1 Heating Condenser
1 DiaPrep Pro (EDTA)

[4007-2001] Duo-Gun Regular Kit
Duo-Gun Handpiece & Charger
GP Obturation Tips (23G, 4/box), Plunger Set, Multi Tool,
Safety Caps (2/pack), Plunger O-Ring Assembly (2/pack),
Scrubbing Brush, Battery, Power Adapter & Cord

[4007-2000] Dia-Duo Regular Kit
Duo-Gun Handpiece, Duo-Pen Handpiece & Duo-Charger,
Heating Condensers .04-50 & .06-50, Disposable Sheath
(200/box), GP Obturation Tips (23G, 4/box), Safety Caps
(2/pack), Plunger Set, Plunger O-Ring Assembly (2/pack),
Multi Tool, Scrubbing Brush, Power Adapter & Cord

[4007-2002] Duo-Pen Regular Kit
Duo-Pen Handpiece & Charger, Heating Condenser
.04-50, & .06-50, Disposable Sheath (200/Box), Power
Adapter & Cord

Duo-Gun™ Features

- Variable temperature settings: 140°C, 170°C, 200°C
- Temperature reaches 200°C within 15 seconds
- Designed for both left-handed and right-handed individuals
- Operate hand trigger with minimal effort
- Strong and durable tips are available in both
23G (#4007-2120) and 25G (#4007-2121)



Duo-Pen™ Features

- Tightly compacts and seals all canals incl. lateral canals
- 3 temperature settings: 170°C, 200°C, and 230°C
- Softens, spreads, cuts, and compacts gutta percha
- Quick heating tip reaches 230°C within 1 second
- 5 color-coded tip sizes: #.04-40, .04-45, .04-50, .06-50, .08-55

Diapex™ Plus

Calcium Hydroxide with Iodoform



BUY 3, GET 1 FREE!*



- 1 Pre-mixed paste eliminates mixing
- 2 Pharmaceutical grade iodoform
- 3 Excellent accessibility to apex
- 4 Antibacterial & bacteriostatic properties
- 5 Multiple indications for use
- 6 Highly radiopaque
- 7 Does not harden, so Ca(OH)₂ remains active

[A1001-501] Diapex Plus Regular Kit

1 x 2g syringe and 20 tips

DIA-ROOT™ Bio Sealer

BioCeramic Root Canal Sealer

BUY 3, GET 1 DIAPEX PLUS (CaOH₂) FREE!



- Excellent biocompatibility
- Highly antibacterial with high pH 12
- Superior bonding strength to dentin & GP
- Zero shrinkage and low solubility (0.95%)
- Washout resistant to prevent micro-leakage
- Optimal flowability of 22mm
- Setting time: 6 hours
- Easy to store with no wasted material

[1003-701] Dia-Root Bio Sealer Regular Kit

1 x 2g syringe, 20 tips, 1 cap