




2023-2024 CATALOG

ENDO TECH is proud to offer **MANI** files.

Advanced products for quality treatment.

25 ¹

The Anatomy of a **MANI** Hand File



- 1** The size and type of the instrument is clearly marked on the handle end.
- 2** Thread eye for floss insertion if required.
- 3** Ergonomically designed handle shaped to comfortably fit to the fingers for ease of use.
- 4** Handle is color coded for quick identification of size.
- 5** Subtle horizontal ribs are incorporated in the handle to prevent slipping when using the instrument.
- 6** Sterilizable, silicon stoppers with directional markers (1mm thick to aid in length adjustments) are on the file. Its shape and color represent the file type and length.
- 7** The highest quality 18-8 stainless steel is treated with a specialized annealing process to produce superior flexibility and hardness.
- 8** The sharpest edge possible for excellent working characteristics and mirror-like finish to prevent corrosion after repeated cleaning and sterilization.
- 9** Taper of 2/100 and 5/100 conforming to the ISO specifications.

Call us today for same-day shipping & personalized service.



ENDO/TECH: OUR STORY	ii
ORDERING, SHIPPING & YOUR SATISFACTION	ii
BURS	1
LAB CUTTING BURS	6
ENDODONTIC ROTARY INSTRUMENTS	7
ENDODONTIC HAND INSTRUMENTS	9
OBTURATION	15
RESTORATIVE MATERIALS	19
BIOCERAMIC MATERIALS	21
SEALERS & TEMP FILLING MATERIALS	22
EDTA AGENTS	23
CALCIUM HYDROXIDE	24
SURGERY	25
ENDODONTIC INSTRUMENTS	26
DENTAL DAM & ACCESSORIES	28
NEEDLES & TIPS	30
EQUIPMENT	33
ULTRASONICS	35
ENDO FOAM & RINGS	42
MATRIX MATERIALS	43
BUR BLOCKS & STANDS	44
TOPICAL AGENTS	45
MIRRORS	46
ACCESSORIES	47



OUR STORY

Endo/Tech was started in 1994 with the intention to give dental professionals a place to get high quality endodontic products at everyday great prices. It was a simple concept and one that has been embraced by the dental community across North America.

The goal of Endo/Tech has always been to treat our customers "right". When you think of it, it is the only way of doing business. Find the "right" products and make sure they are priced "right", all the while providing the "right" service our customers deserve.

When you buy from Endo/Tech you can feel confident that we will not influence your purchasing decisions with high-pressure sales tactics, nor will we inflate our prices to accommodate hidden costs. We are here to help clinicians deliver the highest quality of care by providing and sourcing products as best we can.

Endo/Tech recently celebrated its 29th anniversary thanks to our terrific customers, many of whom have been with us since our inception. Thank you for considering Endo/Tech for your choice for endodontic supplies. We always welcome your comments, questions and suggestions.

HOW TO ORDER

Our office is open weekdays from 7.30am-4.30pm Eastern.

PHONE:	Local Phone: (902) 429-2257	FAX:	Local Fax: (902) 429-5706
	Toll Free Phone: (888) 554-3636		Toll Free Fax: (877) 553-8329

ONLINE:	Website	www.endo-tech.com
	Email	endotech@endo-tech.com

Order online at your convenience, obtain product and manufacturer information and see what's new at Endo/Tech.

SHIPPING

Endo/Tech strives to ensure that all orders are shipped within one business day of the order being placed. Unless otherwise specified by the customer, Endo/Tech ships all orders using a traceable shipping method.

OUR PRODUCT INVENTORY

Endo/Tech maintains a customized inventory based on our customers' needs. As these needs change, we will strive to adjust our inventory to best reflect the product choices of our customers.

ENDO/TECH SATISFACTION GUARANTEE

Endo/Tech aims for 100% customer satisfaction with the products we offer. We also appreciate that needs change and returns do happen. Our return policy is both generous and straightforward.




























Unopened, unused product without an expiry date that was not a special order or discontinued and that has been purchased within the last two years can be returned for credit. Products with an expiry date will be accepted for return within 30 days of purchase. There are no restocking fees unless the manufacturer charges a restocking fee.



























CARBIDE BURS









MANI Standard Operatory Carbides

New and improved standard operatory carbide burs from MANI. Long lasting burs engineered to reduce fracture rates while cutting efficiently.

10 per pack

SHAPE	#	FG
ROUND		
	1/4	MNFG1/4
	1/2	MNFG1/2
	1	MNFG1
	2	MNFG2
	3	MNFG3
	4	MNFG4
	5	MNFG5
	6	MNFG6
	7	MNFG7
	8	MNFG8
INVERTED CONE		
	33 1/2	MNFG331/2
	34	MNFG34
	35	MNFG35
	36	MNFG36
	37	MNFG37
CYLINDER FLAT END PLAIN		
	56	MNFG56
	57	MNFG57
	58	MNFG58
TAPER FLAT END PLAIN		
	169L	MNFG169L
	170	MNFG170
	171	MNFG171
	172	MNFG172
PEAR PLAIN		
	245	MNFG245
	246	MNFG246
	330	MNFG330
	331	MNFG331
	332	MNFG332

SHAPE	#	FG
TAPER FLAT END CROSS CUT		
	556	MNFG556
	557	MNFG557
	558	MNFG558
	559	MNFG559
	560	MNFG560
	699	MNFG699
	700	MNFG700
	701	MNFG701
	701L	MNFG701L
	702	MNFG702
	703	MNFG703
CYLINDER ROUND END PLAIN		
	1156	MNFG1156
	1157	MNFG1157
	1158	MNFG1158
	1159	MNFG1159
CYLINDER ROUND END CROSS CUT		
	1557	MNFG1557
	1558	MNFG1558
	1559	MNFG1559
TRIMMING & FINISHING PLAIN		
	7408	MNFG7408
	7901	MNFG7901
	7902	MNFG7902
	7903	MNFG7903
REMOVAL CROSS CUT		
	1931	MNFG1931
	1957	MNFG1957
	1958	MNFG1958
	1970	MNFG1970














SHAPE	#	FG
SPECIAL SHAPE		
	M16DB	MNFGM16DB
	M14CF	MNFGM14CF
	701X	MNFG701X
	M21DC	MNFGM21DC
	M23DC	MNFGM23DC
	M15EZ	MNFGM15EZ
	CE4CF	MNFGCE4CF
	M15ZB	MNFGM15ZB

CleanCut Sterile Carbide Burs






CleanCut Burs are specifically designed to provide a greater level of blade on tooth contact. Each bur is sterilized and individually packaged.

25 per pack

Friction Grip

SHAPE	#	FG	FG SURG
ROUND			
	1	-	CC1FGSL
	2	CC2FG	CC2FGSL
	4	CC4FG	CC4FGSL
	6	-	CC6FGSL
STRAIGHT/FLAT END CROSS CUT			
	557	CC557FG	CC557FGSL
STRAIGHT/ROUND END CROSS CUT			
	1557	CC1557FG	CC1557FGSL
	1558	CC1558FG	-
STRAIGHT/FLAT END PLAIN CUT			
	56	CC56FG	-
TAPER/FLAT END PLAIN CUT			
	169	CC169FG	-
TAPER/FLAT END CROSS CUT			
	701	CC701FG	-
PEAR			
	245	CC245FG	-
	255SS	CC245SSFG	-
	330	CC330FG	-

Right Angle / Latch

SHAPE	#	RA	RA SURG
ROUND			
	2	CC2RA	CC2RASL
	3	CC3RA	-
	4	CC4RA	CC4RASL
	6	CC6RA	CC6RASL
	8	CC8RA	-

MANI MI Stainless Burs

Specialized stainless steel design for caries removal and conservative dentin removal during endodontic access.

6 per pack

Size	22mm	28mm	34mm	Bur Head Diam.
1/2	MI1/222	MI1/228	-	0.6mm
1	MI1122	MI1128	-	0.8mm
2	MI222	MI228	MI234	1.0mm
3	MI322	MI328	-	1.2mm
4	MI422	MI428	-	1.4mm
5	MI522	MI528	-	1.6mm
6	MI622	MI628	MI634	1.8mm
1-6	MI1-622	MI1-628	-	

MANI Long Tapered Burs

For effective removal of metal posts, fiber cores and canal calcifications. The long taper offers an unobstructed view of the working area.

4 per pack

Size	25mm	28mm
2	MN225	MN228
330	MN33025	MN33028
1557	MN155725	MN155728
MZB	MNMZB25	MNMZB28

1.888.554.3636

2023-2024 PRODUCT CATALOG

SS White Surgical/Endo Length Burs 5 per pack

SHAPE	#	FG	RA
	1	SS1FGSL	-
	2	SS2FGSL	SS2RASL
	4	SS4FGSL	SS4RASL
	6	SS6FGSL	SS6RASL
	8	SS8FGSL	SS8RASL
	557	SS557FGSL	SS557RASL
	558	SS558FGSL	SS558RASL
	1557	SS1557FGSL	SS1557RASL
	1558	SS1558FGSL	-
	700	SS700FGSL	-
	701	SS701FGSL	-
	702	SS702FGSL	SS702RASL
	703	SS703FGSL	SS703RASL

SS White Finishing Burs 5 per pack

	7204	7204 Finish Bur	SS7204FB5
	7404	7404 Finish Bur	SS7404FB5
	7406	7406 Finish Bur	SS7406FB5
	7408	7408 Finish Bur	SS7408FB5
	7642	7642 Finish Bur	SS7642FB5
	7664	7664 Finish Bur	SS7664FB5
	7801	7801 Finish Bur	SS7801FB5
	7802	7802 Finish Bur	SS7802FB5
	7803	7803 Finish Bur	SS7803FB5
	7901	7901 Finish Bur	SS7901FB5
	7902	7902 Finish Bur	SS7902FB5
	CFT3	CFT3 Finish Bur	SSCFT3FB5
	9406	9406 Finish Bur	SSF9406FB

SS White EndoGuide Burs

EndoGuide Burs patented conical-shaped, micro-diameter tip acts as a self-centering guide to create precise endodontic access into canal orifices and around posts. **5 per pack**

Bur #	Head Length (mm)
SSEG1A	2.5
SSEG1	3.5
SSEG2	2.5
SSEG3	1.5
SSEG4	3.5
SSEG5	1.5
SSEG6	2.5
SSEG7	1.5


EndoGuide Molar Kit	SSEGMKIT
EndoGuide Anterior & Bicuspid Kit	SSEGA/BKIT

SS White Standard Operatory Carbides




SHAPE	#	FG 10 Pack	FG 100 Pack	RA 10 Pack
ROUND				
	1/4	SS1/4FG10	SS1/4FG100	-
	1/2	SS1/2FG10	SS1/2FG100	SS1/2RA10
	1	SS1FG10	SS1FG100	SS1RA10
	2	SS2FG10	SS2FG100	SS2RA10
	3	SS3FG10	SS3FG100	SS3RA10
	4	SS4FG10	SS4FG100	SS4RA10
	6	SS6FG10	SS6FG100	SS6RA10
	8	SS8FG10	SS8FG100	SS8RA10
PEAR				
	330	SS330FG10	SS330FG100	SS330RA10
	331	SS331FG10	SS331FG100	SS331RA10
	332	SS332FG10	SS332FG100	-
INVERTED CONE				
	34	SS34FG10	SS34FG100	SS34RA10
	36	SS36FG10	SS36FG100	SS36RA10
	37	SS37FG10	SS37FG100	SS37RA10
	39	SS39FG10	SS39FG100	SS39RA10
STRAIGHT / FLAT END PLAIN				
	56	SS56FG10	SS56FG100	-
	57	SS57FG10	SS57FG100	SS57RA10
	58	SS58FG10	SS58FG100	-
STRAIGHT / FLAT END CROSS CUT				
	556	SS556FG10	SS556FG100	-
	557	SS557FG10	SS557FG100	SS557RA10
	557L	SS557LFG10	SS557LFG100	-
	558	SS558FG10	SS558FG100	SS558RA10
STRAIGHT / ROUND END PLAIN				
	1157	SS1157FG10	SS1157FG100	SS1157RA10
	1158	SS1158FG10	SS1158FG100	-
STRAIGHT / ROUND END CROSS CUT				
	1556	SS1556FG10	SS1556FG100	SS1556RA10
	1557	SS1557FG10	SS1557FG100	SS1557RA10
	1558	SS1558FG10	SS1558FG100	SS1558RA10
TAPER / FLAT END PLAIN				
	169	SS169FG10	SS169FG100	-
TAPER / FLAT END CROSS CUT				
	699	SS699FG10	SS699FG100	-
	700	SS700FG10	SS700FG100	SS700RA10
	701	SS701FG10	SS701FG100	SS701RA10
	701L	SS701LFG10	SS701LFG100	SS701LRA10
	702	SS702FG10	SS702FG100	SS702RA10
TAPER / ROUND END PLAIN				
	1170	SS1170LFG10	SS1170LFG100	SS1170LRA10
AMALGAM PREP				
	245	SS245FG10	SS245FG100	-
	246	SS246FG10	SS246FG100	-
	271	SS271FG10	SS271FG100	-

Visit: www.sswhitedental.com for more information and a full listing of all available bur shapes and sizes.
















SS White Endo Safe End Bur 5 per pack

	ESE014	SS White Endo Safe End Bur	SSESE-014
	ESE018	SS White Endo Safe End Bur	SSESE-018

SS White Great White Surgical Length Burs 5 per pack

GREAT WHITE ROUND			
	GWSL6R	SS White GW 6R Surgical	SSGW6RSL
GREAT WHITE STRAIGHT FISSURE			
	GWSL557	SS White GW 557 Surgical	SSGW557SL
GREAT WHITE ROUND END CROSS CUT			
	GWSL1557	SS White GW 1557 Surgical	SSGW1557SL

SS White Great White Gold Series

SHAPE	#		10 Pack	100 Pack
GREAT WHITE PEAR				
	GW1	SS White GW 1 Bur	SSGW110	SSGW1100
	GW2P	SS White GW 2P Bur	SSGW2P10	-
	GW330	SS White GW 330 Bur	SSGW33010	SSGW330100
GREAT WHITE ROUND				
	GW2R	SS White GW 2R Bur	SSGW2R10	SSGW2R100
	GW4R	SS White GW 4R Bur	SSGW4R10	SSGW4R100
	GW6R	SS White GW 6R Bur	SSGW6R10	SSGW6R100
	GW8R	SS White GW 8R Bur	SSGW8R10	-
GREAT WHITE STRAIGHT FISSURE				
	GW557	SS White GW 557 Bur	SSGW55710	SSGW557100
	GW558	SS White GW 558 Bur	SSGW55810	-
	GW2	SS White GW 2 Bur	SSGW210	SSGW2100
	GW2S	SS White GW 2S Bur	SSGW2S10	SSGW2S10
GREAT WHITE INVERTED CONE				
	GW35	SS White GW 35 Bur	SSGW3510	-
Great White H-Series				
	GWH340	SS White GWH 340 Bur	SSGWH34010	-
	GWH342	SS White GWH 342 Bur	SSGWH34210	-
	GWH342L	SS White GWH 342L Bur	SSGWH342L10	-

SS White Oral Surgery Burs 10 per pack

Shank #1 - 44.5mm

Fits Hall & most standard handpieces

Shank #2 - 51mm

Fits Hall & MicroAire & most standard handpieces

Shank #3 - 65mm

Fits Hall & MicroAire handpieces with long bur guard

Shank #4 - 44.5mm























Fits Stryker 100K & Impaction handpieces

Shank #5 - 59mm

Fits Stryker 100K & Impaction handpieces

Shank #6 - 44.5mm

Fits Stryker 50K & Osteotome Drill J-Notch

		SHANK 1	SHANK 2	SHANK 3	SHANK 4	SHANK 5	SHANK 6
ROUND							
	1	SSO14823	-	-	-	-	-
	2	SSO14826	-	-	-	-	-
	4	SSO14832	-	-	-	-	-
	6	SSO14838	SSO30020	SSO30035	SSO30050	SSO30060	SSO30000
	8	SSO14844	SSO30021	SSO30036	SSO30051	SSO30061	SSO30001
	10	SSO14846	-	-	-	-	-
CROSS CUT FISSURE							
	557	SSO14871	-	-	-	-	-
	558	SSO14874	-	-	-	-	-
	560	SSO14880	-	-	SSO30052	-	SSO30002
	1557	SSO14898	-	-	-	-	-
	1558	SSO14901	-	-	-	-	-
	700	SSO14883	-	-	-	-	-
	701	SSO14886	SSO30022	-	SSO30053	SSO30062	-
	701L	SSO14910	-	-	-	-	-
	702	SSO14889	SSO30023	SSO30037	SSO30054	SSO30063	SSO30003
	702L	SSO14890	SSO30024	-	SSO30055	SSO30064	SSO30004
	1702	SSO14916	SSO30029	-	SSO30048	SSO30069	-
	1702L	SSO14917	SSO30030	-	SSO30049	SSO30070	-
	703	SSO14892	SSO30025	SSO30038	SSO30056	SSO30065	SSO30005
	703L	SSO14893	SSO30026	SSO30039	SSO30057	SSO30066	SSO30006
	1703	SSO14914	SSO30027	SSO30040	SSO30058	SSO30067	SSO30007
	1703L	SSO14915	SSO30028	-	SSO30059	SSO30068	SSO30008

DIAMOND BURS

MANI DiaBurs - Now available in over 400 shapes and grits

Repeat use diamond burs from MANI that provide a superior cutting ability. The electroplating process is strictly controlled by laser measurement inspection and all products exceed ISO standards. For a complete list of the full and complete diamond shapes please visit www.endo-tech.com. **3 to 5 per pack**

MANI coding system reflects shapes of burs. ISO number coding:

TF – Taper Flat End
TR – Taper Round End
TC – Taper Conical End
FO – Flame Ogival End
SF – Straight Flat End

SO – Straight Ogival End
SR – Straight Round End
EX – Special (Extra) Shape
BR – Ball Round Type
BC – Ball Collar Type

DI – Double Inverted Cone
SI – Single Inverted Cone
WR – Wheel Round Edge
CR – CR Inlay Preparation
CD – Children's Dentistry

S – Short Shank
SS – Super Short Shank
C – Coarse
F – Fine
EF – Extra Fine

SC – Super Coarse

YELLOW: Extra Fine (EF) 20-30μ


































RED: Fine (F) 53-63μ





























BLUE: Medium 106-125μ

GREEN: Coarse (C) 125-150μ

BLACK: Super Coarse (SC) 180-212μ

Color Coding - Color Code indication band is provided indicating the size of diamond grit conforming to ISO specifications.

TF SERIES						
	DBTF11			•	•	•
	DBTF12	•	•	•	•	•
	DBTF13			•	•	•
	DBTF14			•	•	
	DBTF20			•	•	
	DBTF21	•	•	•	•	
	DBTF22			•	•	•
	DBTF23			•	•	
	DBTF31			•	•	•
	DBTF41			•	•	
	DBTF42			•	•	
	DBTF43			•	•	
TR SERIES						
	DBTR11	•	•	•	•	
	DBTR12	•	•	•	•	•
	DBTR13	•	•	•	•	•
	DBTR14			•	•	•
	DBTR15			•	•	•
	DBTR19			•	•	•
	DBTR20			•	•	
	DBTR21	•	•	•	•	•
	DBTR24			•	•	
	DBTR25	•	•	•	•	
	DBTR26	•	•	•	•	
TC SERIES						
	DBTC11	•	•	•	•	•
	DBTC21	•	•	•	•	
	DBTC26			•	•	
FO SERIES						
	DBFO11			•	•	
	DBFO21	•	•	•	•	•
	DBFO22	•	•	•	•	•
	DBFO25	•		•	•	
	DBFO27			•	•	•
	DBFO32	•		•	•	•
	DBFO42	•				

SF SERIES						
	DBSF11			•	•	•
	DBSF12			•	•	•
	DBSF13			•	•	•
	DBSF21			•	•	•
	DBSF41			•	•	•
SO SERIES						
	DBSO20			•	•	
	DBSO21			•	•	•
SR SERIES						
	DBSR11			•	•	•
	DBSR12			•	•	•
	DBSR13			•	•	•
EX SERIES						
	DBEX11	•		•	•	•
	DBEX12			•	•	
	DBEX20			•	•	
	DBEX21	•	•	•	•	•
	DBEX24			•	•	
	DBEX26			•	•	•
	DBEX31		•	•	•	
	DBEX41			•	•	
BR SERIES						
	DBBR31			•	•	•
	DBBR41			•	•	•
	DBBR45			•	•	•
	DBBR46			•	•	•
	DBBR49			•	•	
BC SERIES						
	DBBC31			•	•	
	DBBC42			•	•	
	DBBC43			•	•	
DI SERIES						
	DBDI41			•	•	
	DBDI42			•	•	

MANI DiaBurs continued

YELLOW: Extra Fine (EF) 20-30µ


























RED: Fine (F) 53-63µ

















BLUE: Medium 106-125µ

GREEN: Coarse (C) 125-150µ

BLACK: Super Coarse (SC) 180-212µ

Color Coding - Color Code indication band is provided indicating the size of diamond grit conforming to ISO specifications.




SI SERIES						
	DBSI46			•	•	•
	DBSI47			•	•	•
	DBSI48			•	•	•
WR SERIES						
	DBWR13			•	•	•
SHORT SHANK SERIES						
	DBTFS20			•	•	
	DBTFS21			•	•	
	DBTFS22			•		
	DBTFS23			•		
	DBTFS31			•		
	DBTFS41			•		
	DBTRS13			•	•	
	DBTRS21			•	•	
	DBTCS21			•	•	
	DBFOS21			•	•	
	DBSOS20			•	•	
	DBBRS45			•		
	DBBRS46			•	•	
	DBBCS42			•	•	
	DBBCS43			•		
	DBDIS41			•		
	DBSIS46			•		
	DBSIS48			•		
SUPER SHORT SHANK SERIES						
	DBTFS31			•		
	DBTRS21			•	•	
	DBDIS41			•		

CR SERIES	
	DBCR11F
	DBCR11EF
	DBCR12F
	DBCR12EF
	DBCR21F
	DBCR22F
CD SERIES	
	DBCD50F
	DBCD51F
	DBCD52F
	DBCD53F
	DBCD54F
	DBCD55F
	DBCD56F
	DBCD57F
	DBCD58F
	DBCD59F

Please visit www.endo-tech.com for a complete list of all MANI Dia-Burs available





SS White Single Patient Piranha Diamond Burs 25 per pack

*3 most popular burs shown, call for a complete brochure

ROUND			
	801-023C	Piranha Round Coarse	SSP801-023C
FOOTBALL			
	379-023C	Piranha Football Coarse	SSP379-023C
ROUND END TAPER			
	856L-018M	Piranha Round End Taper Medium	SSP856L-018M

SS White Great White Z Diamond Burs 5 per pack

*4 most popular burs shown, call for a complete brochure

ROUND			
	GWZ801-014	GWZ Round	SSGWZ801-014
	GWZ801-018	GWZ Round	SSGWZ801-018
FOOTBALL			
	GWZ379-023	GWZ Football	SSGWZ379-023
ROUND END TAPER			
	GWZ856-018	GWZ Round End Taper	SSGWZ856-018

Clean Cut

LAB CARBIDE BURS

ZIRCONIUM NITRIDE

CleanCut Lab Burs are coated in zirconium nitride to significantly extend the lifespan of the bur. Our coated burs last approximately 10 times longer than standard burs and feature an extremely hard cutting surface of 2,800 Vickers. Other shapes and sizes available as a special order item - please contact us for details.


1 per pack

E SERIES




ITEM #	CCL351.040E	CCL351.060E	CCL79.040E	CCL79.050E	CCL77.023E	CCL77.029E	CCL77.040E	CCL251.060E	CCL73.040E
ISO GRIT	190	190	190	190	190	190	190	190	190
HEAD DIAMETER	040	060	040	050	23	23	23	23	23
HEAD LENGTH	8	13	13	13	5.5	5.5	5.5	5.5	5.5
REC. SPEED	30K	20K	30K	30K	40K	40K	40K	40K	40K
BRASSER #	H351E.11.040	H351E.11.060	H79E.11.040	H79E.11.050	H77E.11.023	H77E.11.029	H77E.11.040	H73.11.040	H77E.11.023
EDENTA #	6210.040HP	5210.060HP	5710.040HP	-	7710.023HP	7710.029HP	-	-	7710.023HP

GSQ SERIES



ITEM #	CCL79.070GSQ	CCL251.060GSQ	CCL261.023GSQ	CCL351.060GSQ
ISO GRIT	172	172	172	172
HEAD DIAMETER	070	060	023	060
HEAD LENGTH	14	14	11.5	13
REC. SPEED	20K	20K	40K	20K
BRASSER #	H79GSQ.11.070	H251GSQ.11.060	H261GSQ.11.023	H351GSQ.11.060
EDENTA #	-	7275.060HP	0175.023HP	5275.060HP

EF SERIES



ITEM #	CCL79.040EF	CCL251.060EF
ISO GRIT	140	140
HEAD DIAMETER	040	060
HEAD LENGTH	13	14
REC. SPEED	30K	20K
BRASSER #	H79EF.11.040	H251EF.11.060
EDENTA #	5720.040HP	5420.060HP

NICKEL TITANIUM FILES



CleanCut Reciprocating File System

The CleanCut Reciprocating File System is designed to safely and efficiently shape the root canal space by using a reciprocating motion. Available in 4 sizes and 3 different lengths. For use with existing reciprocating handpieces. Equivalent sizing and usage to Wave One Gold files. Made in the USA.

3 per pack / assorted 4 pack

	21mm	25mm	31mm
SMALL	CCRS21	CCRS25	CCRS31
PRIMARY	CCRP21	CCRP25	CCRP31
MEDIUM	CCRM21	CCRM25	CCRM31
LARGE	CCRL21	CCRL25	CCRL31
15-45 ASSORTED	-	CCRA25	-



CleanCut Tapered File System

The CleanCut Tapered File System is available in the same sizes, lengths and tapers as the Dentsply® ProTaper Gold® file system and Edge Endo® Edge Taper™ file system. Made in the USA. **6 per pack**

	19mm	21mm	25mm	31mm
SX	CCTSX	-	-	-
S1	-	CCTS121	CCTS125	CCTS131
S2	-	CCTS221	CCTS225	CCTS231
F1	-	CCTF121	CCTF125	CCTF131
F2	-	CCTF221	CCTF225	CCTF231
F3	-	CCTF321	CCTF325	CCTF331
F4	-	CCTF421	CCTF425	CCTF431
F5	-	CCTF521	CCTF525	CCTF531
SX-F3 Asst	-	CCTA21	CCTA25	CCTA31

JIZAI File System by **MANI**



Each JIZAI file is fully heat treated, highly flexible and designed to cut smoothly and efficiently. The cross sectional shape and flute pitch helps to prevent over instrumentation and file binding. Each file can be pre bent and will hold its shape. Made in Japan. **3 per pack**

	18mm	21mm		25mm		31mm	
Size		.04	.06	.04	.06	.04	.06
25.14 ORIFICE OPENER	MNJ251418	-	-	-	-	-	-
25	-	MNJ250421	MNJ250621	MNJ250425	MNJ250625	MNJ250431	MNJ250631
30	-	MNJ300421	MNJ300621	MNJ300425	MNJ300625	MNJ300431	MNJ300631
35	-	MNJ350421	MNJ350621	MNJ350425	MNJ350625	MNJ350431	MNJ350631
40	-	MNJ400421	-	MNJ400425	-	MNJ400431	-
45	-	MNJ450421	-	MNJ450425	-	MNJ450431	-
50	-	MNJ500421	-	MNJ500425	-	MNJ500431	-
25.04 / 25.06 / 35.04 Kit	-	MNJST21		MNJST25		MNJST31	

JIZAI Glider Glide Path File by **MANI**

The JIZAI Glider is ideal for establishing a glide path and easing the transition to greater tapered NiTi shaping files. JIZAI Gliders are ideally suited to fit into the MANI JIZAI Endodontic System but can be used as part of a practitioners existing filing protocol. JIZAI Gliders are manufactured from an optimized R-Phase/Austenitic NiTi alloy to allow the file to exhibit excellent cutting ability despite the finer tip size. Made in Japan. **3 per pack**

	21mm	25mm	31mm
Size	.04	.04	.04
13	MNJ130421	MNJ130425	MNJ130431

Gutta Percha Remover Files by **MANI**

MANI GPR Files offer an innovative approach to removing gutta percha without the need for repurposing NiTi rotary files or using excessive amounts of solvent. **4 per pack**

#	Size	Taper	Length	Material
1S	#70	0.04	16mm	Stainless
2S	#50	0.04	18mm	Stainless
3N	#40	0.04	21mm	NiTi
4N	#30	0.04	21mm	NiTi
Asst	-	0.04	-	x2 SS x2 NiTi

NICKEL TITANIUM FILES

Pac-File
(Twisted / Vortex Style)
Files

With the proprietary treatment of the NiTi alloy, Pac-File™ files do not have the spring back effect as regular NiTi files. These files stay more centered in the canal and conform better to the canal morphology reducing the risk of straightening curved canals, ledge formation, transportation, and perforations.

6 per pack

	21mm		25mm		31mm	
Size	.04	.06	.04	.06	.04	.06
15	PDF150421	PDF150621	PDF150425	PDF150625	PDF150431	PDF150631
20	PDF200421	PDF200621	PDF200425	PDF200625	PDF200431	PDF200631
25	PDF250421	PDF250621	PDF250425	PDF250625	PDF250431	PDF250631
30	PDF300421	PDF300621	PDF300425	PDF300625	PDF300431	PDF300631
35	PDF350421	PDF350621	PDF350425	PDF350625	PDF350431	PDF350631
40	PDF400421	PDF400621	PDF400425	PDF400625	PDF400431	PDF400631
45	PDF450421	PDF450621	PDF450425	PDF450625	PDF450431	PDF450631
50	PDF500421	PDF500621	PDF500425	PDF500625	PDF500431	PDF500631
55	PDF550421	-	PDF550425	-	PDF550431	-
60	PDF600421	-	PDF600425	-	PDF600431	-
15-40	PDFA15-400421	PDFA15-400621	PDFA15-400425	PDFA15-400625	PDFA15-400431	PDFA15-400631

One File G
(Reciprocating) Files

The OneFile G system is designed to safely and efficiently shape canals in a reciprocating motion, and is compatible with the same motor, handpiece and setting as WaveOne® Gold Files.

3 per pack (5 per assorted pack)

	21mm	25mm	31mm
15.02 GLIDE PATH	PDO1521	PDO1525	PDO1531
20.07 SMALL	PDO2021	PDO2025	PDO2031
25.07 PRIMARY	PDO2521	PDO2525	PDO2531
35.06 MEDIUM	PDO3521	PDO3525	PDO3531
45.05 LARGE	PDO4521	PDO4525	PDO4531
15-45 ASSORTED	-	PDOA25	-

DRILLS, REAMERS & CARRIERS

Gates Drills

Gates Drills are designed with a flame-shaped, boring head with a non-cutting guide tip on a slender shaft. **6 per pack**



MANI

Size	28mm	32mm	38mm	diameter mm
1	GD128	GD132	GD138	0.50
2	GD228	GD232	GD238	0.70
3	GD328	GD332	GD338	0.90
4	GD428	GD432	GD438	1.10
5	GD528	GD532	GD538	1.30
6	GD628	GD632	GD638	1.50
1-6	GD1-628	GD1-632	GD1-638	

Peeso Reamers

Peeso Reamers are designed with long, sharp flutes which provide rapid widening of the coronal portion of the canal. Often used for post preparation. **6 per pack**



MANI

Size	28mm	32mm	38mm	diameter
1	PR128	PR132	PR138	0.70
2	PR228	PR232	PR238	0.90
3	PR328	PR332	PR338	1.10
4	PR428	PR432	PR438	1.30
5	PR528	PR532	PR538	1.50
6	PR628	PR632	PR638	1.70
1-6	PR1-628	PR1-632	PR1-638	

Paste Carriers / Lentulo Spirals

Manufactured from spring stainless steel with tightly wound safety coils at the base of the shank. Tapered design to improve material distribution and air bubble release during use.

4 per pack



MANI

Size / # of Grooves	21mm	25mm
● 25 (1)	PC2521	PC2525
● 30 (2)	PC3021	PC3025
● 35 (3)	PC3521	PC3525
● 40 (4)	PC4021	PC4025
25-40	PC25-4021	PC25-4025

MANI K Files / H Files / Reamers



These high quality files are manufactured from high quality stainless steel blanks. **6 per pack**

MANI K Files

	Stainless Steel with rubber stops		
Size	21mm	25mm	31mm
06	0621K	0625K	0631K
08	0821K	0825K	0831K
10	1021K	1025K	1031K
12	1221K	1225K	-
15	1521K	1525K	1531K
17	1721K	1725K	-
20	2021K	2025K	2031K
22	2221K	2225K	-
25	2521K	2525K	2531K
27	2721K	2725K	-
30	3021K	3025K	3031K
32	3221K	3225K	-
35	3521K	3525K	3531K
37	3721K	3725K	-
40	4021K	4025K	4031K
45	4521K	4525K	4531K
50	5021K	5025K	5031K
55	5521K	5525K	5531K
60	6021K	6025K	6031K
70	7021K	7025K	7031K
80	8021K	8025K	8031K
90	9021K	9025K	9031K
100	10021K	10025K	10031K
110	11021K	11025K	11031K
120	12021K	12025K	12031K
130	13021K	13025K	13031K
140	14021K	14025K	14031K
15-40	A15-4021K	A15-4025K	A15-4031K
45-80	A45-8021K	A45-8025K	A45-8031K
90-140	A90-14021K	A90-14025K	A90-14031K

MANI H Files

	Stainless Steel with rubber stops		
Size	21mm	25mm	31mm
06	-	-	-
08	0821H	0825H	0831H
10	1021H	1025H	1031H
12	1221H	1225H	-
15	1521H	1525H	1531H
17	1721H	1725H	-
20	2021H	2025H	2031H
22	2221H	2225H	-
25	2521H	2525H	2531H
27	2721H	2725H	-
30	3021H	3025H	3031H
32	3221H	3225H	-
35	3521H	3525H	3531H
37	3721H	3725H	-
40	4021H	4025H	4031H
45	4521H	4525H	4531H
50	5021H	5025H	5031H
55	5521H	5525H	5531H
60	6021H	6025H	6031H
70	7021H	7025H	7031H
80	8021H	8025H	8031H
90	9021H	9025H	9031H
100	10021H	10025H	10031H
110	11021H	11025H	11031H
120	12021H	12025H	12031H
130	13021H	13025H	13031H
140	14021H	14025H	14031H
15-40	A15-4021H	A15-4025H	A15-4031H
45-80	A45-8021H	A45-8025H	A45-8031H
90-140	A90-14021H	A90-14025H	A90-14031H

MANI Reamers

	Stainless Steel with rubber stops		
Size	21mm	25mm	31mm
06	0621R	0625R	0631R
08	0821R	0825R	0831R
10	1021R	1025R	1031R
12	1221R	1225R	-
15	1521R	1525R	1531R
17	1721R	1725R	-
20	2021R	2025R	2031R
22	2221R	2225R	-
25	2521R	2525R	2531R
27	2721R	2725R	-
30	3021R	3025R	3031R
32	3221R	3225R	-
35	3521R	3525R	3531R
37	3721R	3725R	-
40	4021R	4025R	4031R
45	4521R	4525R	4531R
50	5021R	5025R	5031R
55	5521R	5525R	5531R
60	6021R	6025R	6031R
70	7021R	7025R	7031R
80	8021R	8025R	8031R
90	9021R	9025R	9031R
100	10021R	10025R	10031R
110	11021R	11025R	11031R
120	12021R	12025R	12031R
130	13021R	13025R	13031R
140	14021R	14025R	14031R
15-40	A15-4021R	A15-4025R	A15-4031R
45-80	A45-8021R	A45-8025R	A45-8031R
90-140	A90-14021R	A90-14025R	A90-14031R

Premier K Files / H Files / Reamers

6 per pack



PREMIER K FILES

	Stainless Steel with a silicone stop		
Size	21mm	25mm	31mm
06	PM0621K	PM0625K	PM0631K
08	PM0821K	PM0825K	PM0831K
10	PM1021K	PM1025K	PM1031K
12	-	-	-
15	PM1521K	PM1525K	PM1531K
17	-	-	-
20	PM2021K	PM2025K	PM2031K
22	-	-	-
25	PM2521K	PM2525K	PM2531K
27	-	-	-
30	PM3021K	PM3025K	PM3031K
32	-	-	-
35	PM3521K	PM3525K	PM3531K
37	-	-	-
40	PM4021K	PM4025K	PM4031K
45	PM4521K	PM4525K	PM4531K
50	PM5021K	PM5025K	PM5031K
55	PM5521K	PM5525K	PM5531K
60	PM6021K	PM6025K	PM6031K
70	PM7021K	PM7025K	PM7031K
80	PM8021K	PM8025K	PM8031K
90	PM9021K	PM9025K	PM9031K
100	PM10021K	PM10025K	PM10031K
110	-	-	-
120	-	-	-
130	-	-	-
140	-	-	-
15-40	PMA15-4021K	PMA15-4025K	PMA15-4031K
45-80	PMA45-8021K	PMA45-8025K	PMA45-8031K
90-140	-	-	-

PREMIER H FILES

	Stainless Steel with a silicone stop		
Size	21mm	25mm	31mm
08	PM0821H	PM0825H	PM0831H
10	PM1021H	PM1025H	PM1031H
12	-	-	-
15	PM1521H	PM1525H	PM1531H
17	-	-	-
20	PM2021H	PM2025H	PM2031H
22	-	-	-
25	PM2521H	PM2525H	PM2531H
27	-	-	-
30	PM3021H	PM3025H	PM3031H
32	-	-	-
35	PM3521H	PM3525H	PM3531H
37	-	-	-
40	PM4021H	PM4025H	PM4031H
45	PM4521H	PM4525H	PM4531H
50	PM5021H	PM5025H	PM5031H
55	PM5521H	PM5525H	PM5531H
60	PM6021H	PM6025H	PM6031H
70	PM7021H	PM7025H	PM7031H
80	PM8021H	PM8025H	PM8031H
90	PM9021H	PM9025H	PM9031H
100	PM10021H	PM10025H	PM10031H
110	-	-	-
120	-	-	-
130	-	-	-
140	-	-	-
15-40	PMA15-4021H	PMA15-4025H	PMA15-4031H
45-80	PMA45-8021H	PMA45-8025H	PMA45-8031H
90-110	-	-	-
90-140	-	-	-

PREMIER REAMERS

	Stainless Steel with a silicone stop		
Size	21mm	25mm	31mm
06	PM0621R	PM0625R	PM0631R
08	PM0821R	PM0825R	PM0831R
10	PM1021R	PM1025R	PM1031R
12	-	-	-
15	PM1521R	PM1525R	PM1531R
17	-	-	-
20	PM2021R	PM2025R	PM2031R
22	-	-	-
25	PM2521R	PM2525R	PM2531R
27	-	-	-
30	PM3021R	PM3025R	PM3031R
32	-	-	-
35	PM3521R	PM3525R	PM3531R
37	-	-	-
40	PM4021R	PM4025R	PM4031R
45	PM4521R	PM4525R	PM4531R
50	PM5021R	PM5025R	PM5031R
55	PM5521R	PM5525R	PM5531R
60	PM6021R	PM6025R	PM6031R
70	PM7021R	PM7025R	PM7031R
80	PM8021R	PM8025R	PM8031R
90	PM9021R	PM9025R	PM9031R
100	PM10021R	PM10025R	PM10031R
110	-	-	PM11031R
120	-	-	PM12031R
130	-	-	PM13031R
140	-	-	PM14031R
15-40	PMA15-4021R	PMA15-4025R	PMA15-4031R
45-80	PMA45-8021R	PMA45-8025R	PMA45-8031R
90-140	-	-	-

*allow 2-3 weeks for shipping

*allow 2-3 weeks for shipping

Micro Files

Micro Files from MANI are ideal for finding and maintaining a visual field under a microscope. They are particularly useful when working on posterior and other regions where instruments are difficult to manipulate.

K Micro Files are ideal for finding canals and initial instrumentation of the root canal space. H Micro Files are for enlarging a root canal, removing gutta percha, debris and pastes. Flare Micro Files feature a 0.05 tapered design for enlarging a root canal and checking the flare preparation.

Micro Files are used with the MANI Endo Holder. Perfectly balanced and featuring an ergonomic textured design that comes to a round end for easy attachment of an apex locator file clip, the Endo Holder is a useful addition to your endodontic armamentarium.

K Type Micro Files:	0.02 Taper	23mm Length	Sizes #08-40	6 per pack
H Type Micro Files:	0.02/0.04 Taper	23mm Length	Sizes #08-40	6 per pack
F Type Micro Files:	0.05 Taper	23mm Length	Sizes #15-30	6 per pack



MANI Micro Files



Stainless Steel			
Size	K Files	H Files	Flare Files
08	MNM0802K	MNM0804H	-
10	MNM1002K	MNM1004H	-
15	MNM1502K	MNM1502H	MNM1505F
20	MNM2002K	MNM2002H	MNM2005F
25	MNM2502K	MNM2502H	MNM2505F
30	MNM3002K	MNM3002H	MNM3005F
35	MNM3502K	MNM3502H	-
40	MNM4002K	MNM4002H	-
15-40	MNMA15-4002K	MNMA15-4002H	-

FOR CALCIFIED CANALS

MANI D-Finders

The D-shaped cross section allows unprecedented initial access and navigation of calcified and curved canals.

6 per pack



Size	21mm	25mm
● 08	0821D	0825D
● 10	1021D	1025D
● 12	1221D	1225D
○ 15	1521D	1525D

MANI Glide Finders

Glide Finders are designed for patency and creating a glide path. Compared to other files designed for calcified canals and path finding, Glide Finders have higher resistance to buckling, flexibility and cutting efficiency.

6 per pack

Size	21mm	25mm
● 08	0821G	0825G
● 10	1021G	1025G
● 12	1221G	1225G
○ 15	1521G	1525G

NiTi HAND FILES

DiaDent NiTi Hand K Files

Reamers and Hedstrom also available.

6 per pack

Size	21mm	25mm	31mm
● 08	DNT0821K	DNT0825K	DNT0831K
● 10	DNT1021K	DNT1025K	DNT1031K
○ 15	DNT1521K	DNT1525K	DNT1531K
● 20	DNT2021K	DNT2025K	DNT2031K
● 25	DNT2521K	DNT2525K	DNT2531K
● 30	DNT3021K	DNT3025K	DNT3031K
● 35	DNT3521K	DNT3525K	DNT3531K
● 40	DNT4021K	DNT4025K	DNT4031K
15-40	DNT15-4021K	DNT15-4025K	DNT15-4031K
45-80	DNT45-8021K	DNT45-8025K	DNT45-8031K

SPECIALTY HAND FILES

MANI Flare Files

Flare Files offer good flexibility and cutting ability and are very effective for curved root canals. Flare Files reduce the number of instruments needed as compared to the conventional step-back preparation. Triangle cross-section with a .05 taper. **6 per pack**



Size	21mm	25mm
○ 15	1521FF	1525FF
● 20	2021FF	2025FF
● 25	2521FF	2525FF
● 30	3021FF	3025FF
● 35	3521FF	3525FF
● 40	4021FF	4025FF
○ 45	4521FF	4525FF
● 50	5021FF	5025FF
● 55	5521FF	5525FF
● 60	6021FF	6025FF
15-40	A15-4021FF	A15-4025FF

MANI Sec-O Files (Safe End K and H)

The rounded tip of the Sec-O File slides on the outer wall of a curved canal, ensuring safe preparation. The safety tip of these files is designed to prevent ledging and perforation during root canal preparation.

6 per pack

**K MANI Sec-O K**

Size	21mm	25mm
● 08	0821SK	0825SK
● 10	1021SK	1025SK
○ 15	1521SK	1525SK
● 20	2021SK	2025SK
● 25	2521SK	2525SK
● 30	3021SK	3025SK
● 35	3521SK	3525SK
● 40	4021SK	4025SK
○ 45	4521SK	4525SK
● 50	5021SK	5025SK
● 55	5521SK	5525SK
● 60	6021SK	6025SK
● 70	7021SK	7025SK
● 80	8021SK	8025SK
15-40	A15-4021SK	A15-4025SK
45-80	A45-8021SK	A45-8025SK

**H MANI Sec-O H**

Size	21mm	25mm
● 08	-	-
● 10	1021SH	1025SH
○ 15	1521SH	1525SH
● 20	2021SH	2025SH
● 25	2521SH	2525SH
● 30	3021SH	3025SH
● 35	3521SH	3525SH
● 40	4021SH	4025SH
○ 45	4521SH	4525SH
● 50	5021SH	5025SH
● 55	5521SH	5525SH
● 60	6021SH	6025SH
● 70	7021SH	7025SH
● 80	8021SH	8025SH
15-40	A15-4021SH	A15-4025SH
45-80	A45-8021SH	A45-8025SH

FLEXIBLE FILES

MANI RT Files

RT Files have a unique cross sectional design for greater flexibility and efficient cutting ability. It provides higher performance in a super flexible, stainless steel instrument. **6 per pack**



Size	21mm	25mm	31mm
○ 15	1521RT	1525RT	1531RT
● 20	2021RT	2025RT	2031RT
● 25	2521RT	2525RT	3131RT
● 30	3021RT	3025RT	3031RT
● 35	3521RT	3525RT	3531RT
● 40	4021RT	4025RT	4031RT
○ 45	4521RT	4525RT	4531RT
● 50	5021RT	5025RT	5031RT
● 55	5521RT	5525RT	5531RT
● 60	6021RT	6025RT	6031RT
● 70	7021RT	7025RT	7031RT
● 80	8021RT	8025RT	8031RT
15-40	A15-4021RT	A15-4025RT	A15-4031RT
45-80	A45-8021RT	A45-8025RT	A45-8031RT

MANI Flexile Files

The triangular shape of the Flexile File improves the flexibility and cutting efficiency, especially when preparing a curved root canal. **6 per pack**



Size	21mm	25mm
○ 15	1521F	1525F
● 20	2021F	2025F
● 25	2521F	2525F
● 30	3021F	3025F
● 35	3521F	3525F
● 40	4021F	4025F
15-40	A15-4021F	A15-4025F

MANI Steri-Endo & Bur Guards

The Steri-Endo Guard enables easy removal of debris or dirt from both hand files and rotary files. Fully autoclavable when loaded with files, only the sponges should be considered disposable. 16 files per guard. **1 per pack**

The MANI 12 hole and 22 hole Bur Guards are fully autoclavable and fit both Latch and Friction Grip burs. The height of each Bur Guard can be adjusted to fit both standard length and surgical length burs. **1 per pack**



Blue Steri-Endo Guard	MNEGB
Yellow Steri-Endo Guard	MNEGY
White Steri-Endo Guard	MNEGW

Blue 12 Hole Bur Guard	MNBG12B
Yellow 12 Hole Bur Guard	MNBG12Y
White 12 Hole Bur Guard	MNBG12W

Blue 22 Hole Bur Guard	MNBG22B
Yellow 22 Hole Bur Guard	MNBG22Y
Green 22 Hole Bur Guard	MNBG22G

Finger Spreaders

Finger Spreaders have a conical tip and are tapered to match accessory sized gutta percha points. Used for lateral condensation techniques. **6 per pack**

MANI

	Size	21mm	25mm
○	15 / XF	FS1521	FS1525
●	20 / FF	FS2021	FS2025
●	25 / MF	FS2521	FS2525
●	30 / F	FS3021	FS3025
●	35 / FM	FS3521	FS3525
●	40 / M	FS4021	FS4025
	15-40	FS15-4021	FS15-4025



Premier

	Size	21mm	25mm
	15 / XF	PM1521FS	PM1525FS
	20 / FF	PM2021FS	PM2025FS
	25 / MF	PM2521FS	PM2525FS
	30 / F	PM3021FS	PM3025FS
	35 / FM	PM3521FS	PM3525FS
	40 / M	PM4021FS	PM4025FS
	15-40	PMA15-4021FS	PMA15-4025FS

Finger Pluggers

Finger Pluggers have a flat tip and a .02 taper, used for lateral and vertical condensation. **6 per pack**

MANI

	Size	21mm	25mm
○	15	FP1521	FP1525
●	20	FP2021	FP2025
●	25	FP2521	FP2525
●	30	FP3021	FP3025
●	35	FP3521	FP3525
●	40	FP4021	FP4025
	15-40	FP15-4021	FP15-4025



Premier

	Size	21mm	25mm
	15	PM1521FP	PM1525FP
	20	PM2021FP	PM2025FP
	25	PM2521FP	PM2525FP
	30	PM3021FP	PM3025FP
	35	PM3521FP	PM3525FP
	40	PM4021FP	PM4025FP
	15-40	PMA15-4021FP	PMA15-4025FP

Barbed Broaches

MANI - metal handles - **6 per pack**



MANI

	Size	
○	XXXF	BB121
●	XXF	BB221
●	XF	BB321
●	F	BB421
●	M	BB521
●	C	BB621
	Assorted	BBA2-421

GUTTA PERCHA

ISO Standard Sized Gutta Percha

Laser inspected, latex-free and hand rolled. Spill-proof package is plastic with 6 compartments.

Spillproof Pack **120 per pack**

Endo/Tech RGS

	Size	Spillproof Color Coded
●	25	GP25
●	30	GP30
●	35	GP35
●	40	GP40
○	45	GP45
●	50	GP50
●	55	GP55
●	60	GP60
●	70	GP70
●	80	GP80


DiaDent

	Size	Spillproof Color Coded
●	8	GPD08
●	10	GPD10
○	15	GPD15
●	20	GPD20
●	25	GPD25
●	30	GPD30
●	35	GPD35
●	40	GPD40
○	45	GPD45
●	50	GPD50
●	55	GPD55
●	60	GPD60
●	70	GPD70
●	80	GPD80
	15-40	GPD415-40
	45-80	GPD45-80

Also available in mm-marked

Gutta Percha Obturator Pellets

100 per jar



GPEO



GPDO

Accessory Gutta Percha

Spillproof Pack or Feather Tip **100 per pack**

Endo/Tech RGS

	Size	Spillproof Color Coded	Feather Tip Color Coded
●	FF	GPFF	GPFFFT
●	MF	GPMF	GPMFFT
●	F	GPFF	GPFFFT
●	FM	GPFM	GPFMFT
●	M	GPM	GPMFT
○	ML	GPML	-
●	L	GPL	-

DiaDent

	Size	Spillproof Color Coded	Feather Tip Color Coded
○	XF	GPDXF	GPDXFFT
●	FF	GPdff	GPdffFT
●	MF	GPDmf	GPDmfFT
●	F	GPdf	GPdffFT
●	FM	GPDfm	GPDfmFT
●	M	GPDm	GPDmFT
○	ML	GPDml	GPDmlFT
●	L	GPDl	GPDlFT
●	XL	GPDXL	GPDXLFT
	XF-M	GPDXF-M	GPDXF-MFT

Also available in mm-marked

Extra Long Gutta Percha

Oversized Dia-VET gutta percha points are available to order. They are designed for unusually long canals and are available in 35mm, 45mm and 60mm lengths with a standard 0.02 taper.

Greater Taper Gutta Percha ISO GT

Complement: Profile ISO Rotary Files and K3/Twisted Files. Spillproof box. **60 per pack**

Endo/Tech RGS

Size	.04 Taper	.06 Taper
15	GPEI1504	GPEI1506
20	GPEI2004	GPEI2006
25	GPEI2504	GPEI2506
30	GPEI3004	GPEI3006
35	GPEI3504	GPEI3506
40	GPEI4004	GPEI4006
45	GPEI4504	GPEI4506
50	GPEI5004	GPEI5006
55	GPEI5504	GPEI5506
60	GPEI6004	GPEI6006
70	GPEI7004	GPEI7006
80	GPEI8004	GPEI8006
15-40	GPEIA15-4004	GPEIA15-4006

DiaDent

Size	.04 Taper	.06 Taper
10	GPD11004	GPD11006
15	GPD11504	GPD11506
20	GPD12004	GPD12006
25	GPD12504	GPD12506
30	GPD13004	GPD13006
35	GPD13504	GPD13506
40	GPD14004	GPD14006
45	GPD14504	GPD14506
50	GPD15004	GPD15006
55	GPD15504	GPD15506
60	GPD16004	GPD16006
70	GPD17004	GPD17006
80	GPD18004	GPD18006
15-40	GPDIA15-4004	GPDIA15-4006
45-80	GPDIA45-8004	GPDIA45-8006

Also available in mm-marked

Dia-PRO W Gutta Percha

Complement: Wave One Files. **60 per pack**

Small	GPDPWS
Primary	GPDPWP
Large	GPDPWL
Assorted	GPDPWA

Dia-PRO T NEXT Gutta Percha

Complement: ProTaper NEXT Files
Millimeter Marked, length - 28mm.
60 per pack

X2	GPDPX2
X3	GPDPX3
X4	GPDPX4
X5	GPDPX5
X2-X5	GPDPX2-5
X2-X3	GPDPX2-3

Dia-PRO T Gutta Percha

Complement: ProTaper Files
Millimeter Marked, length - 28mm.
60 per pack

F1	GPDPF1
F2	GPDPF2
F3	GPDPF3
F4	GPDPF4
F5	GPDPF5
F1/F2/F3	GPDPF1-3
F4/F5	GPDPF4-5
F1-F5	GPDPFALL

DiaDent Dia-PRO W GOLD Gutta Percha

Complement: Wave One Gold Files. **60 per pack**

Small	GPDPWGS
Primary	GPDPWGP
Medium	GPDPWGM
Large	GPDPWGL
Assorted	GPDPWGA

DiaDent Gutta Percha PRO

Complement: Profile Series 29 Files. **60 per pack**

Size	.04 Taper	.06 Taper
#2	GPDPRO042	GPDPRO062
#3	GPDPRO043	GPDPRO063
#4	GPDPRO044	GPDPRO064
#5	GPDPRO045	GPDPRO065
#6	GPDPRO046	GPDPRO066
#7	GPDPRO047	GPDPRO067
#8	GPDPRO048	GPDPRO068
#9	GPDPRO049	GPDPRO069
#10	GPDPRO0410	GPDPRO0610
#2-7	GPDPRO04A2-7	GPDPRO06A2-7

DiaDent Dia-SYSTEM GTP Gutta Percha

Complement: ProSystem GT, ProFile GT Files, GT Series X, EdgeFile & EdgeEvolve. **60 per pack**

Size	.04 Taper	.06 Taper	.08 Taper	.10 Taper	.12 Taper
20	GPDSG2004	GPDSG2006	GPDSG2008	GPDSG2010	GPDSG2012
25	GPDSG2504	GPDSG2506	GPDSG2508	GPDSG2510	-
30	GPDSG3004	GPDSG3006	GPDSG3008	GPDSG3010	-
35	GPDSG3504	GPDSG3506	GPDSG3508	GPDSG3510	-
40	GPDSG4004	GPDSG4006	GPDSG4008	GPDSG4010	-

PAPER POINTS

ISO Standard Sized Paper Points

Made from specially formulated, highly absorbent papers and hand rolled to ensure superior absorbency. Cell package points remain sterile until opened. Spill-proof package is plastic with 6 compartments.

Spillproof Pack **200 per pack**

Endo/Tech RGS

	Size	Spillproof Color Coded
○	15	PP15
●	20	PP20
●	25	PP25
●	30	PP30
●	35	PP35
●	40	PP40
○	45	PP45
●	50	PP50
●	55	PP55
●	60	PP60
●	70	PP70
●	80	PP80

DiaDent

	Size	Spillproof Color Coded	Cells Color Coded
●	8	PPD08	PPD08C
●	10	PPD10	PPD10C
○	15	PPD15	PPD15C
●	20	PPD20	PPD20C
●	25	PPD25	PPD25C
●	30	PPD30	PPD30C
●	35	PPD35	PPD35C
●	40	PPD40	PPD40C
○	45	PPD45	PPD45C
●	50	PPD50	PPD50C
●	55	PPD55	PPD55C
●	60	PPD60	PPD60C
●	70	PPD70	PPD70C
●	80	PPD80	PPD80C
	15-40	PPDA15-40	PPDA15-40C
	45-80	PPDA45-80	PPDA45-80C

Also available in mm-marked

Extra Long Paper Points

Oversized Dia-VET paper points are available to order. They are designed for unusually long canals and are available in 35mm, 45mm and 60mm lengths with a standard 0.02 taper.

Accessory Paper Points

Spillproof Pack or Cell Pack. **200 per pack**


Endo/Tech RGS

	Size	Spillproof Non Color Coded
●	XF	PPXF
●	F	PPF
●	M	PPM
●	C	PPC
●	XC	PPXC

DiaDent

	Size	Spillproof Color Coded	Cells Color Coded	Feather Tip Color Coded
○	XXF	PPDXXF	PPDXXFC	PPDXXFFT
●	XF	PPDXF	PPDXFC	PPDXFFT
●	F	PPDF	PPDFC	PPDFFT
●	M	PPDM	PPDMC	PPDMFT
●	C	PPDC	PPDCC	PPDCFT
●	XC	PPDXC	PPDXCC	PPDXCFT
	XXF-XC	PPDAXXF-XC	PPDAXXF-XCC	PPDAXXF-XCFT

Also available in mm-marked



Greater Taper Paper Points ISO GT

Spillproof Pack. **100 per pack**



Endo/Tech RGS



Size	.04 Taper	.06 Taper
15	PPEI1504	PPEI1506
20	PPEI2004	PPEI2006
25	PPEI2504	PPEI2506
30	PPEI3004	PPEI3006
35	PPEI3504	PPEI3506
40	PPEI4004	PPEI4006
45	PPEI4504	PPEI4506
50	PPEI5004	PPEI5006
55	PPEI5504	PPEI5506
60	PPEI6004	PPEI6006
70	PPEI7004	PPEI7006
80	PPEI8004	PPEI8006

DiaDent

Size	.04 Taper	.06 Taper
10	PPDI1004	PPDI1006
15	PPDI1504	PPDI1506
20	PPDI2004	PPDI2006
25	PPDI2504	PPDI2506
30	PPDI3004	PPDI3006
35	PPDI3504	PPDI3506
40	PPDI4004	PPDI4006
45	PPDI4504	PPDI4506
50	PPDI5004	PPDI5006
55	PPDI5504	PPDI5506
60	PPDI6004	PPDI6006
70	PPDI7004	PPDI7006
80	PPDI8004	PPDI8006
15-40	PPDIA15-4004	PPDIA15-4006
45-80	PPDIA45-8004	PPDIA45-8006

Also available in mm-marked

Dia-PRO W Paper Points

Complement: Wave One Files.
100 per pack

Small	PPDPWS
Primary	PPDPWP
Large	PPDPWL
Assorted	PPDPWA

Dia-PRO T NEXT Paper Points

Complement: ProTaper NEXT Files
100 per pack

X2	PPDPX2
X3	PPDPX3
X4	PPDPX4
X5	PPDPX5
X2-X5	PPDPX2-5
X2-X3	PPDPX2-3

Dia-PRO T Paper Points

Complement: ProTaper Files.
100 per pack

F1	PPDPF1
F2	PPDPF2
F3	PPDPF3
F4	PPDPF4
F5	PPDPF5
F1/F2/F3	PPDPF1-3
F4/F5	PPDPF4-5
F1-F5	PPDPFALL

DiaDent Dia-PRO W GOLD Paper Points

Complement: Wave One Gold Files.
100 per pack

Small	PPDPWGS
Primary	PPDPWGP
Medium	PPDPWGM
Large	PPDPWGL
Assorted	PPDPWGA

Dia-PRO Paper Points

Complement: Profile Series 29 Files.
100 per pack



Size	.04 Taper	.06 Taper
#2	PPDPRO042	PPDPRO062
#3	PPDPRO043	PPDPRO063
#4	PPDPRO044	PPDPRO064
#5	PPDPRO045	PPDPRO065
#6	PPDPRO046	PPDPRO066
#7	PPDPRO047	PPDPRO067
#8	PPDPRO048	PPDPRO068
#9	PPDPRO049	PPDPRO069
#10	PPDPRO0410	PPDPRO0610
#2-7	PPDPRO04A2-7	PPDPRO06A2-7

DiaDent Dia-SYSTEM GTP Paper Points

Complement: ProSystem GT, ProFile GT Files, GT Series X, EdgeFile & EdgeEvolve.
100 per pack

Size	.04 Taper	.06 Taper	.08 Taper	.10 Taper	.12 Taper
20	PPDSG2004	PPDSG2006	PPDSG2008	PPDSG2010	PPDSG2012
25	PPDSG2504	PPDSG2506	PPDSG2508	PPDSG2510	-
30	PPDSG3004	PPDSG3006	PPDSG3008	PPDSG3010	-
35	PPDSG3504	PPDSG3506	PPDSG3508	PPDSG3510	-
40	PPDSG4004	PPDSG4006	PPDSG4008	PPDSG4010	-

MANI RESTORATIVE MATERIALS

NOW IN
CANADA

MANI RESTORATIVE PRODUCTS

Now available in Canada. MANI RESTORATIVE composites are the perfect combination of Japanese Standards and German Technology.

MANI Micro - MicroHybrid Composite

MANI MICRO is a universal light-curing filling material with a high filler content that's easy to use. MANI MICRO offers excellent esthetics and strength for all your dental composite needs.



The special combination of macro and micro fillers gives this material its outstanding physical properties and makes it easy to polish. Thanks to these excellent qualities, this material is extremely plaque resistant. Scientifically proven to be especially biocompatible and gentle to soft tissue.

MANI MICRO, Available in A1, A2, A3, A3.5, A4, B1 & B2 & CHAMELEON

MRM

MANI NANO - NanoFill Composite

MANI NANO - True nano technology

MANI NANO has an ultra fine filler for exceptionally homogenous restorations. Featuring excellent mechanical properties for optimized working time and results and available in a wide range of shades with high physical values.



MANI NANO, Available in A1-4, B1-3, C2 & C3

MRN

MANI BULK - BulkFill Composite

MANI BULK is a posterior restorative material, visible, light curing optimized to create posterior restorations more simply and faster. MANI BULK provides excellent strength and low wear for durability with up to a 4 mm depth of cure. MANI BULK is both flowable and highly viscous.

Available in a universal chameleon shade for all your needs.



MANI BULK, Available in Colour Matching CHAMELEON

MRBUS

ACTIVA BioActive



ACTIVA BioACTIVE Restorative

ACTIVA BioACTIVE-RESTORATIVE is a highly esthetic, bioactive composite that delivers all the advantages of glass ionomers in a strong, resilient, resin matrix that will not chip or crumble. It chemically bonds to teeth, seals against microleakage and releases more calcium, phosphate and fluoride. Activa is more bioactive than glass ionomers, and is more durable and fracture resistant than composites.

BioACTIVE Restorative Starter Kit, Available in A1, A2, A3	PUVR
BioACTIVE Restorative Single Refill, Available in A1, A2, A3, A3.5	PUVR1
BioACTIVE Restorative Value Refill, Available in A1, A2, A3, A3.5	PUVR2



ACTIVA BioACTIVE Base Liner

ACTIVA BioACTIVE-BASE/LINER delivers all the advantages of glass ionomers in a strong, resilient resin matrix that will not chip or crumble. It chemically bonds to teeth, seals against microleakage and releases calcium, phosphate and fluoride. It is more bioactive than glass ionomers, and is more durable and fracture resistant than composites, glass ionomers and RMGIs.

BioACTIVE Base Liner Single Pack	PUVB1
BioACTIVE Base Liner Value Pack	PUVB2



ACTIVA BioACTIVE Cement

ACTIVA BioACTIVE-CEMENT stimulates mineral apatite formation and the natural remineralization process at the material-tooth interface that knits together the restoration and the tooth, penetrates and fills micro-gaps, reduces sensitivity, guards against secondary caries, and seals margins against microleakage and failure.

BioACTIVE Cement Single Pack, Available in Translucent or A2 Opaque	PUVC1
BioACTIVE Cement Value Pack, Available in Translucent or A2 Opaque	PUVC2



ACTIVA Kids

ACTIVA KIDS is an esthetic, BioACTIVE resin restorative that stimulates apatite formation and the natural remineralization process with release and recharge of calcium, phosphate and fluoride. It is an opaque white shade ideally suited for pediatric dentistry.

Kids Starter Kit, Available in Opaque	PUVKP
Kids Single Refill, Available in Opaque	PUVK1P
Kids Value Refill, Available in Opaque	PUVK2P

BIOCERAMIC MATERIALS



Avalon BioMed Bioceramic NeoSEALER Flo

NeoSEALER Flo is a bioactive bioceramic root canal sealer with superior handling properties, promoting hydroxyapatite formation to trigger the healing process. Unlike conventional sealers, NeoSEALER Flo is biocompatible, antimicrobial, dimensionally stable and RESIN-FREE.

NeoSEALER Refill: 2.2g Syringe (No Tips)	AVSRS
NeoSEALER Starter Kit: 2.2g + 20 Flex Flo Tips	AVSSK
NeoSEALER Professional Kit: 4.4g (2 x 2.2gm) + 40 Flex Flo Tips	AVSPK
NeoSEALER Refill Bundle: 10 x 2.2g (22gm, No Tips)	AVSRSB
NeoSEALER Professional Bundle: 5 x 4.4g + 200 Flex Flo Tips (22gm)	AVSPKB
NeoSEALER Flex Flo Tips 20Pk	AVFT20



Avalon BioMed Bioceramic NeoPUTTY

NeoPUTTY™ is a bioactive bioceramic premixed root and pulp treatment with superior handling properties, promoting hydroxyapatite formation to trigger the healing process.

NeoPUTTY 1.2g Professional Kit (Approx 16 doses)	AVP
NeoPUTTY BUNDLE 5x 1.2g (Approx 80 doses)	AVPB



Avalon BioMed Bioceramic NeoMTA 2

Stainproof, tricalcium silicate-based bioactive cement that can be used universally for vital pulp and other endodontic indications in primary and permanent teeth. NeoMTA 2 bioceramic triggers the healing process.

NeoMTA 2 x3 2.5g Smart Bundle	AVN2B
NeoMTA 2 2.5g Kit	AVN2



DiaDent DiaRoot Bio Sealer

DiaRoot Bio Sealer is a pre-mixed bioceramic calcium silicate-based MTA sealer. Its non-shrinking feature and strong bonding ability make it the perfect sealer for permanent obturation and suitable for all obturation techniques with gutta percha.

Bio Sealer Regular Kit 2g Syringe + 20 Tips	D1003-701
---	-----------



NovaCal Pulp Capping Material with MTA Fillers

NovaCal releases calcium-ions which promotes the formation of tertiary dentine. The high pH level has bactericidal effects, supports healing and protects against hypersensitivity. Due to its thixotropic behaviour, NovaCal can be applied very precisely, even in deep cavity preps.

NovaCal 1g Syringe + 10 Tips	NOC1
------------------------------	------

SEALERS



DiaDent Dia-ProSeal

Dia-ProSeal is a 2 component epoxy resin based system similar in design to AH Plus. High pH value due to the addition of calcium hydroxide. Excellent radiopacity, setting time is 3-4 hours.

Regular Kit Contents: 1 16g syringe, mixing pad, spatula.

Intro Kit Contents: 1 4g syringe, mixing pad, spatula.

Dia-ProSeal Regular Kit	D1003-201
Dia-ProSeal Intro Kit	D1003-202



Kerr/Sybron Pulp Canal Sealer EWT

Kerr Pulp Canal Sealer is a zinc oxide eugenol sealer that has a powder base and liquid catalyst. It is a non-toxic and non-irritating formulation that is very radiopaque and has a long history of clinical success. EWT features a working time of greater than 6 hours on the pad.

Contents: 10.5g powder, 4ml liquid, 2 dispensing scoops.

Kerr/Sybron Pulp Canal Sealer EWT	KSEWT
-----------------------------------	-------

TEMPORARY FILLING MATERIAL



DiaDent e-Temp

e-Temp adheres to dentin, expands slightly during setting and hardens rapidly to ensure a hermetical seal. Moisture-activated, sets within 25 minutes or less.

e-Temp 30g Jar	D2004-1001
----------------	------------



3M ESPE Cavit

Cavit is a self-curing, radiopaque temporary restorative. Offering easy and complete sealing of endodontic access preparations. 28g per jar.

Cavit-P (Pink - Hard Set)	CV-PJ
Cavit-W (White - Medium Set)	CV-WJ
Cavit-G (Gray - Soft Set)	CV-PGJ



DiaDent DiaTemp LC Temporary Restorative Material

DiaTemp is a high-viscosity light-cured restorative material for temporary and intermediate restorations for inlays and onlays. It can also be used to reline prefabricated crowns and bridges or to seal implant screw channels.

Flowable Light-Cured Temporary Restorative Material (Yellow)	DA2003-2201
Flowable Light-Cured Temporary Restorative Material (Blue)	DA2003-2101
Light-Cured Temporary Restorative Material (Yellow)	DA2003-1201
Light-Cured Temporary Restorative Material (Blue)	DA2003-1101

EDTA



DiaDent Dia-Prep Pro

Ethylenediamine Tetraacetic Acid (EDTA) 17%.

Urea Peroxide 3.5%.

Polypropylene Glycol & Polyethylene Glycol.

Contents: 2 x 6g Dia-Prep Plus syringes and 10 disposable tips.

DiaDent Dia-Prep Pro	D2002-4100
----------------------	------------



Pulpdent EDTA 17% Solution

This 17% solution of EDTA has been buffered to a neutral pH.

EDTA 120ml bottle	PUEDA120
EDTA 480ml bottle	PUEDA480



Pulpdent Prep-Rite

17% viscous EDTA gel, contains peroxide and lubricant.

Kit Contents: 4 x 5g syringes.

Prep-Rite RC Kit	PUPR
------------------	------



Premier RC-Prep

RC-Prep's formulation of glycol, urea peroxide and EDTA in a water-soluble base helps remove calcification and lubricates the canal for efficient instrumentation.

Available in the following forms:	
RC-Prep 18g Jar	PMRCP
RC-Prep 227g Jar	PMRCP227

Calcium Hydroxide



DiaDent Diapex Plus

Diapex is used in root canal filling therapy in either temporary or permanent procedures. The addition of Iodoform to the calcium hydroxide improves radiopacity and adds an antibacterial agent to the paste.

Kit Contents: 2g of pre-mixed paste in a pre-filled syringe + 20 disposable tips.

DiaDent Diapex Complete Kit	D1001-501
-----------------------------	-----------



Pulpdent Multi-Cal

All Purpose, Non-Setting Calcium Hydroxide Paste. Bactericidal - Bacteriostatic - pH >12 - Radiopaque. Multi-Cal is a smooth, creamy calcium hydroxide preparation recommended for all clinical applications where calcium hydroxide is indicated.

Kit Contents: 4 x 1.2ml syringes + 8 applicator tips.

Refill Syringe: 3ml syringe.

Refill Tips: Black, 22ga x 1/2", prebent tips, 100 per pack.

Multi-Cal Kit	PUMCK
Multi-Cal Refill Syringe	PUMC
Multi-Cal Refill Tips Black / Dark Blue	PUMCT B / DB



Pulpdent Paste

Pulpdent Paste is the original premixed calcium hydroxide methylcellulose pulpal dressing with barium sulfate. It is indicated for direct pulp capping, pulpotomy and as an intracanal dressing between office visits.

Kit Contents: 3ml syringe + 24 needles (18 gauge blunt, sterile).

Pulpdent Paste Kit	PUPPK
--------------------	-------



Pulpdent TempCanal Enhanced

Flows through 27ga x 25mm endodontic irrigation needles with 2-side-vents for controlled placement and uniform coating of canal walls, and with a closed end to prevent overfills. Easy flow formula, non drying paste, pH>12, radiopaque.

Kit Contents: 4 x 1.2ml syringes + 20 needles (27 gauge, 2 side vents, sterile).

3ml Syringe: 1 x 3ml luer lock syringe, no needles.

TempCanal Enhanced 4 Syringe Kit	PUTE4
TempCanal Enhanced 3ml Syringe	PUTE3



J&J Ethicon Sutures

Call for other sizes or visit woundclosure.ca.

1654G 4-0 Chromic Gut, 18" in Length, P-3 Needle	JJ1654G
687G 5-0 Chromic Gut, 18" in Length, P-3 Needle	JJ687G
1637G 4-0 Chromic Gut, 18" in Length, PS-2 Needle	JJ1637G
1636G 5-0 Chromic Gut, 18" in Length, PS-3 Needle	JJ1636G
1644G 4-0 Plain Gut, 18" in Length, P-3 Needle	JJ1644G
612G 3-0 Plain Gut, 18" in Length, X-1 Needle	JJ612G
1676G 5-0 Black Silk, 18" in Length, PS-2 Needle	JJ1676G
1677G 4-0 Black Silk, 18" in Length, PS-2 Needle	JJ1677G
J495G 5-0 Coated Vicryl, 18" in Length, PS-2 Needle	JJJ495G
J496G 4-0 Coated Vicryl, 18" in Length, PS-2 Needle	JJJ496G
J499G 6-0 Coated Vicryl, 18" in Length, PS-3 Needle	JJJ499G
Y464G 4-0 Monocryl Dyed, 18" in Length, P-3 Needle	JJY464G
Y463G 5-0 Monocryl Dyed, 18" in Length, P-3 Needle	JJY463G
Y513G 4-0 Monocryl Dyed, 18" in Length, PS-2 Needle	JJY513G

G = 12/pk H = 36/pk

Please allow a 1 week delivery time for non-stocked Ethicon sutures. For sale in Canada only.



Miltex Sterile Surgical Blades

15C angled blade from Miltex. Available in both carbon and stainless steel. Carbon blades are sharper but dull more quickly. **100 per box**

Miltex 15C Stainless Steel Sterile Surgical Blades, 100 per pack	MX4-315C
--	----------



Paragon Sterile Surgical Blades

Precision cutting edge design. Convenient sterile package. Available in Stainless Steel and Carbon Steel, in the following sizes:

10 - 11 - 12 - 15. **100 per box**

Paragon Disposable Sterile Scalpels

Convenient single use and comfort fit handle design. Precision molding of blade into handle with a replaceable blade guard for complete protection. Available in the following sizes:

10 - 11 - 12 - 15. **10 per box**



Hemostatic Products

Spongostan is a sterile, water-insoluble, malleable, porcine gelatin absorbable sponge intended for hemostatic use by applying to a bleeding surface.

J&J Spongostan - Absorbable Gelatin Sponge Hemostat, 24 per pack	JJSP
J&J Surgifoam 2cm x 6cm x 7mm, 12 per pack	JJ1972
J&J Surgicel 0.5 inches x 2 inches, 12 per pack	JJ1955

Endodontic Instruments

EXPLORERS

Endo/Tech RGS Explorers

Crafted in Germany from the highest quality surgical stainless steel. Features a comfortable round serrated solid handle.



Endo/Tech RGS Explorer, DG16	EIDG16
Endo/Tech RGS Explorer, Fine, DG16F	EIDG16F

Standard tipped DG16 are angled 70° and 55° from handle. Ultra fine tips have a 45° and 50° angle.

Endo/Tech Black Chromium Plated RGS Explorers

The special chromium plating significantly hardens the stainless steel to improve sharpness and strength while maintaining tactile feedback.



Endo/Tech RGS Black Explorer, DG16	EIBDG16
Endo/Tech RGS Black Explorer, Fine, DG16F	EIBDG16F



Premier DG16 RC Explorer DE	PMEDG16
-----------------------------	---------



Premier DG23/6 RC Explorer DE	PMEDG23/6
-------------------------------	-----------



Premier Finder 2	PMF2
------------------	------

SPOON EXCAVATORS



Endo/Tech RGS Excavator 31L DE 1.0	EIE31L
------------------------------------	--------



Endo/Tech RGS Excavator 32L DE 1.5	EIE32L
------------------------------------	--------



Endo/Tech RGS Excavator 33L DE 2.0	EIE33L
------------------------------------	--------



Premier Excavator 31L DE	PME31L
--------------------------	--------



Premier Excavator 31LS DE	PME31LS
---------------------------	---------



Premier Excavator 32L DE	PME32L
--------------------------	--------



Premier Excavator 33L DE	PME33L
--------------------------	--------

Endo/Tech carries the complete line of Kohler, Miltex and Premier Endo instruments. Please call for pricing.

PLUGGERS



DiaDent 35/70 Dia-Kondensor	D551-101
DiaDent 40/80 Dia-Kondensor	D551-102
DiaDent 50/100 Dia-Kondensor	D551-103
DiaDent 60/120 Dia-Kondensor	D551-104
DiaDent Dia-Kondensor Complete Set	D551-191



Premier NC Plugger 1/3	PMPN1/3
------------------------	---------



Premier NC Plugger 5/7	PMPN5/7
------------------------	---------



Premier NC Plugger 9/11	PMPN9/11
-------------------------	----------

SPREADERS



Premier D11 Spreader SE	PMSD11
-------------------------	--------



Premier D11T Spreader SE	PMSD11T
--------------------------	---------



Premier NiTi D11 Spreader SE	PMSND11
------------------------------	---------



Premier NiTi D11T Spreader SE	PMSND11T
-------------------------------	----------

CARRIERS / PLIERS



Endo/Tech RGS Locking Plier, Dual Grooved	EIL
---	-----

Endo/Tech RGS Non-Locking Plier, Dual Grooved	EINL69
---	--------

Endo/Tech RGS Non-Locking Plier, Serrated	EINLS
---	-------



Miltex Self Locking Serrated Plier	MX66-652
------------------------------------	----------



Miltex Self Locking Grooved Plier	MX66-656
-----------------------------------	----------



Premier Endo Non-Locking Plier	PMNLEP
--------------------------------	--------

OTHER



Endo/Tech RGS Glick #1	EIGL1
------------------------	-------



Endo/Tech RGS Perio Probe	EIP23W
---------------------------	--------



Premier Root Canal Glick #1	PMG1
-----------------------------	------



Premier Root Canal Glick #2	PMG2
-----------------------------	------

DENTAL DAMS & ACCESSORIES



TOP
SELLER



isolate Latex Dental Dam

isolate Dental Dam is manufactured with a rolled sheeting process similar to extrusion. This process provides superior quality control over the more commonly used dipping process.

isolate Dental Dam is minimally powdered and has a unique non-reflective glare surface. Blue colored latex is low protein and powder free.

5X5 **52 per pack** 6X6 **36 per pack**

5 X 5

	HEAVY	MEDIUM	THIN
Blue	i5HB	-	-
Green	-	i5MG	i5TG
Lite Sand	-	i5MLS	-
Purple	-	i5MP	-
Rainbow	-	-	-

6 X 6

	HEAVY	MEDIUM	THIN
Blue	i6HB	i6MB	-
Green	i6HG	i6MG	i6TG
Lite Sand	i6HLS	i6MLS	-
Purple	-	i6MP	-

isolate Latex-Free Dental Dam

isolate Latex-Free Dental Dam is manufactured from synthetic polyisoprene and is completely latex free. isolate Latex-Free Dental Dam provides superior elasticity and strength.

5X5 **30 per pack** 6X6 **20 per pack**
5X5 Bulk **100 per pack** 6X6 Bulk **75 per pack**

5 X 5

	MEDIUM
Blue	i5NLMB
Blue Bulk	i5BULK

6 X 6

	MEDIUM
Blue	i6NLMB
Blue Bulk	i6BULK

Puldent Kool-Dam

Heatless Liquid Dam & Block Out Resin. Apply Kool-Dam on gingiva or tooth surfaces prior to power bleaching, sandblasting or other procedures requiring intraoral protection, or use to block out undercuts prior to taking impressions.

Kit Contents: 2 x 3ml syringes and 20 applicator tips

Bulk Kit Contents: 10 x 3ml syringes, no applicator tips

Kool-Dam Kit	PUKDK
Kool-Dam Bulk Kit	PUKDB

Dental Dam Frames 1 per pack

Miltex Metal Frame	MX7625
Young Plastic Frame	YDFP

DENTAL DAM ACCESSORIES

 Premolars
 Anteriors
 Molars

isolate Dental Dam Forceps



isolate Forcep - Ivory

iGFI



isolate Forcep - Stoke

iGFS

isolate Dental Dam Punches



isolate Punch - Ainsworth

iGPA



isolate Punch - Ivory

iGPI

isolate Dental Dam Clamps

Featuring improved flexibility, superior strength and accurate spring tension.



Clamp 00

i00



Clamp 0

i0



Clamp 1

i1



Clamp 1T

i1T



Clamp 2

i2



Clamp 2A

i2A



Clamp 2AT

i2AT



Clamp 3

i3



Clamp 4

i4



Clamp 7

i7



Clamp 8

i8



Clamp 8A

i8A



Clamp W8A

iW8A



Clamp 9

i9



Clamp 9T

i9T



Clamp 12A

i12A



Clamp 13A

i13A



Clamp 14

i14



Clamp 14A

i14A



Clamp 14T

i14T



Clamp 26

i26



Clamp 56

i56



Clamp 209

i209



Clamp 212

i212



Painless Steel Injection Needles

Painless Steel Dental Injection Needles are disposable, sterile needles for the delivery of anesthetic. Each needle has a three-edge lancet grinding and a siliconized canula for smooth entry and comfort to patient. Each hub has a bevel marking and is threaded to ensure easy and safe handling and to prevent leakage.

Benefits:

Ultra sharp and siliconized needle slides smoothly into tissue.

No more worrying about the direction of deflection because of the bevel mark.

Threaded hub fits safely and securely on the syringe for easy handling.

Precise base gluing ensures a secure and accurate fit of the canula.

100 tips per pack

25Ga Short Injection Needles (Red)	TRPS25S
25Ga Long Injection Needles (Red)	TRPS25L
27Ga Short Injection Needles (Yellow)	TRPS27S
27Ga Long Injection Needles (Yellow)	TRPS27L
30Ga Extra Short Injection Needles (Blue)	TRPS30XS
30Ga Short Injection Needles (Blue)	TRPS30S



Luer-Lock Endo Irrigation Syringes

Ideal for site irrigation, compatible with all luer style tips and hubs.

Disposable and non-sterile. Contains no latex components. **50-100 per pack**

1.2cc	100pk	PDIS1.2
3cc	100pk	PDIS3
6cc	50pk	PDIS6
12cc	100pk	PDIS12



Curved Utility Syringe


Disposable, non-sterile, ideal for irrigation or intra-oral impression material delivery. **50 per pack**

Curved Utility Syringe - 12 cc	PDUS
--------------------------------	------

Luer Lock Handpiece

Combined with 39mm extra long suction tips and a low volume saliva ejector, the Luer Lock Handpiece provides an enhanced field of view and precise suction in areas difficult to access. **1 per pack**

Luer Lock Handpiece	TRFC3100
---------------------	----------



Sterile Blunt End Tips

Sterile, bendable tips with a rounded cannula, compatible with all luer style syringes. **100 tips per pack**

23mm Long Blunt Tip, 27Ga, Blue	TR6030
23mm Long Blunt Tip, 30Ga, Green	TR6031



39mm Extra Long Suction Tips

These 39mm luer lock tips are ideal for direct suction in areas that are difficult to access, e.g. during conventional or surgical endodontic interventions.

100 tips per pack

14Ga 39mm Round End Suction Green Tip	TR6480
16Ga 39mm Round End Suction White Tip	TR6481
18Ga 39mm Round End Suction Pink Tip	TR6482
20Ga 39mm Round End Suction Yellow Tip	TR6483



PD125

PD1357

PD1169

Delivery Tips

Delivery tips are excellent choices for the precise and controlled application of dental materials. **100 tips per pack**

Delivery Tips - Brush Pad	PD125
Delivery Tips - Curved With Brush	PD1169
Delivery Tips - Micro Brush	PD1357



Liner

Flow

Flow Sealant

Sealant Cement

Etch

Pre-Bent Tips

Disposable, blunt end pre-bent needle tips for use with etchants, resins and flowable composites. **100 tips per pack**

Pre-Bent Blunt Tip, 18Ga, Pink	PD227
Pre-Bent Blunt Tip, 20Ga, Yellow	PD228
Pre-Bent Blunt Tip, 20Ga, Black	PD229
Pre-Bent Blunt Tip, 22Ga, Gray	PD230
Pre-Bent Blunt Tip, 25Ga, Blue	PD231




Precision Molded Tips

Precision molded plastic tips ideal for irrigation, suction and delivery. Luer lock connection.

20 tips per pack

Pre-Bent Tips, Clear, 20pk	D1001-320
----------------------------	-----------



Irrigation Tips

Bendable 1" monoject style (notched) tips, color coded for easy identification. **100 tips per pack**


Pac-Dent, 23Ga	PD238
Pac-Dent, 27Ga	PD239
Pac-Dent, 30Ga	PD232



Closed End Side Port Irrigation Tips

Designed with a closed end cannula and a side window dispersal point for effective irrigation. **100 tips per pack**

OptiProbe, 23Ga	PDITO23
OptiProbe, 27Ga	PDITO27
OptiProbe, 30Ga	PDITO30



Pre-Bent Closed End Side Port Irrigation Tips

Designed with a closed end cannula and a side window dispersal point for effective irrigation. The pre-bent design is adjustable to facilitate ease of access. **50 tips per pack**

	Single Side Port	Double Side Port
27Ga, 21mm	PD217	PD215
31Ga, 21mm	PD219	PD221
27Ga, 27mm	PD218	PD216
31Ga, 27mm	PD220	PD222



Double Side Vented Irrigation Tips

Featuring two lateral side vents and a closed end cannula for effective irrigation. Each tip features thin wall construction for enhanced flow rates. Each tip is individually sterilized. **25-40 tips per pack**

Irrigation Tip, 23Ga, Blue, 40pk	TR6490
Irrigation Tip, 25Ga, Orange, 40pk	TR6491
Irrigation Tip, 27Ga, Grey, 25pk	TR6402
Irrigation Tip, 30Ga, Yellow, 25pk	TR6403
Irrigation Tip, 31Ga, Purple, 25pk	TR6415



3/4" Evacu-N-Irrigation Tips

3/4", blunt-end, bendable, compatible with all luer style syringes.

100 tips per pack

Evacu-N-Irrigation Blunt Tip, 18Ga, Pink	PD233
Evacu-N-Irrigation Blunt Tip, 19Ga, Grey	PD234
Evacu-N-Irrigation Blunt Tip, 20Ga, Yellow	PD235
Evacu-N-Irrigation Blunt Tip, 22Ga, Black	PD236
Evacu-N-Irrigation Blunt Tip, 25Ga, Blue	PD237

Dia-Duo Obturation System

Dia-Duo is the new advanced cordless obturation system from DiaDent that consists of both Duo-Gun BackFill Obturation System and Duo-Pen Warm Vertical Compaction System. Contact us for full technical specifications.



Dia-DUO Complete Obturation System	D4007-2000
Duo-Pen Vertical Compaction Device	D4007-2002
Duo-Gun Backfill Obturation Device	D4007-2001
Duo-Pen 40.04 Compaction Tip	D4007-2070
Duo-Pen 45.04 Compaction Tip	D4007-2071
Duo-Pen 50.04 Compaction Tip	D4007-2072
Duo-Pen 50.06 Compaction Tip	D4007-2073
Duo-Pen 55.08 Compaction Tip	D4007-2074
Duo-Gun 23Ga Needle Tips	D4007-2120
Duo-Gun 25Ga Needle Tips	D4007-2121

Obtura MaxPack System

The Obtura Spartan MaxPack System can be used to remove excess gutta percha when performing lateral condensation techniques, create an endodontic post space, or compact gutta percha to prepare for warm gutta percha injection techniques. Contact us for full technical specifications.



Obtura III Max Heated Gutta Percha System	OS824-700
MaxPack Plugger 35.04	OS824-711
MaxPack Plugger 40.04	OS824-712
MaxPack Plugger 45.04	OS824-713
MaxPack Plugger 55.06	OS824-714
MaxPack Plugger 55.08	OS824-715
MaxPack Plugger 55.10	OS824-716
MaxPack Plugger 60.12	OS824-717

Obtura III Max Heated Gutta Percha System

Getting a comprehensive fill in primary canals is a good thing. But when you take that final x-ray and you see gutta percha has also flowed into lateral and accessory canals, fins, loops and isthmuses, you know you've achieved total obturation excellence. Those are the results you can expect with the Obtura III Max, the instrument that heats gutta percha to its ideal viscosity for MAXimum obturation success. Backed by over 25 years of trusted clinical performance. Contact us for full technical specifications.



Obtura III Max Heated Gutta Percha System	OS823-800
Obtura III Gutta Percha Pellets	OS822-602
Obtura III Delivery Needles 20Ga (5Pk)	OS823-620
Obtura III Delivery Needles 23Ga (5Pk)	OS823-623
Obtura III Delivery Needles 25Ga (5Pk)	OS823-635



EvoFill Duo Cordless Obturation System

With an innovative electric motor that prevents hand fatigue, EvoFill Duo offers fast heating and controlled gutta percha extrusion with precise 3D fills and reliable results. EvoFill has a user replaceable battery. Includes x1 Gutta Percha Cartridge and x1 heating condenser tip.

EvoFill Duo Cordless Obturation System	D4005-1001
EvoFill Gutta Percha Cartridge (20 per box)	D4005-2010
Extra Fine Condenser Tip	D4005-2070
Fine Condenser Tip	D4005-2071
Fine Medium Condenser Tip	D4005-2072
Medium Condenser Tip	D4005-2073
Medium Large Condenser Tip	D4005-2074



Newtron Booster



Newtron Booster B. LED



Newtron Booster B. LED XS

Acteon Newtron Ultrasonic Units

The Newtron range of ultrasonic generators from Acteon are equipped with exclusive and patented Newtron technology, which brings more preservation, efficacy and comfort to any procedure.

Preservation

- Controlled Vibrations - More precise treatment and gentle on teeth
- Tip Quality - Anatomy and enamel preservation
- Total Irrigation Control - No heating, safer for tissue

Efficacy

- Frequency Adjustment - Maximum efficiency no matter the load applied to the unit
- Power Regulation - Consistent performance, no pressure needed
- Powerful Cavitation - Fragmentation and elimination of deposits, cleaning and disinfecting effect

Comfort

- Linear & Regular Vibrations - Painless treatment for the patient, less fatigue for the clinician and preserves the sense of tactile feedback.

Newtron Booster	ACF60206
Newtron Booster B. LED	ACF61106
Newtron Booster B. LED XS	ACF62106

Canal Access Preparation



CAP1 (1 per pack)

ACF88181

Micro-blade tip length 12mm, taper 6%

The CAP1 tip is used with the active lateral part for:

- Finishing walls and polishing.
- Removing temporary cement and dentinal residues.
- Removing dentin overhangs.

The CAP1 tip has a non-active end to prevent the risk of perforating the pulp chamber floor.



CAP2 (1 per pack)

ACF88182

Micro-blade tip, length 9mm, taper 5%

The CAP2 tip has active lateral part and extremity and is used by sweeping method to remove dentin bridges.

- Location of the MB2 (2nd mesiobuccal canal) and search for hidden canals.
- Preparation of the pulp chamber.
- Removal of the dentin layer which may hide the access to the MB2 canal.

The micro-blades are less aggressive than diamond and their coating makes these tips very durable.



CAP3 (1 per pack)


ACF88183

Micro-blade tip, length 8mm taper 6%

The CAP3 tip has a very pointed extremity indicated for:

- Locating and opening the calcified canals.
- Fragmenting calcifications or pulp stones in the pulp chamber.
- Loosening fiber posts.
- Locating accessory canals.

Due to its very sharp point, the CAP3 tip must be handled with care (visual aids recommended).



ET18D (1 per pack)

ACF88017

Diamond-coated steel tip 76µm, length 18mm, taper 5%

The ET18D tip is a diamond-coated tip for:

- Finishing the access cavity.
- Removing dentin overhangs, calcifications and filling materials.



ETBD (1 per pack)

ACF88020

Diamond-coated ball tip, length 20mm, taper 5%

Diamond-coated ball tip for searching for canals and locating calcified canals.



Canal Irrigation



IRRI 20/21mm Irrigating Files (4 per pack)	ACF43807
IRRI 20/25mm Irrigating Files (4 per pack)	ACF43808
IRRI 25/21mm Irrigating Files (4 per pack)	ACF43805
IRRI 25/25mm Irrigating Files (4 per pack)	ACF43806


Passive ultrasonic irrigation (PUI) files

Irrisafe™ files are used for irrigation once the root canal has been prepared. The file is inserted 2mm short of the working length and activated by performing withdrawal movements to flush the debris and the smear layer upwards.

Irrisafe™ safely eliminates the smear layer, dentine debris and bacteria from the root canal. Its blunt tip prevents any risk of perforating the apex or the canal walls.



K10/21mm (4 per pack)	ACF43710
K10/25mm (4 per pack)	ACF43712
K15/21mm (4 per pack)	ACF43715
K15/25mm (4 per pack)	ACF43717
K25/21mm (4 per pack)	ACF43725
K25/25mm (4 per pack)	ACF43727
K30/21mm (4 per pack)	ACF43730
K30/25mm (4 per pack)	ACF43732



Files of different lengths and diameters, taper 2%

K files can be used for irrigation, removing calcified dentine and gutta percha, and retrieving broken instruments.

Obturation



SO4 (1 per pack)	ACF88009
-------------------------	----------

Fine condenser, length 40mm, taper 4%

The SO4 tip is designed for lateral condensation of gutta percha by heating effect. It is used dry, without irrigation.

Canal Retreatment




ET20 (1 per pack)

ACF88011

Retreatment tip, length 20mm, taper 6%

The ET20 steel tip is used in the 1st coronal third:

- Extraction of filling material, silver points, broken instruments.
- Removal of debris and the smear layer.



ET20D (1 per pack)

ACF88013

Diamond-coated retreatment tip, 30 μ m, length 20mm, taper 5%

The ET20D tip in diamond-coated steel is used in the 1st coronal third to remove very hard materials by brushing the walls.

The diamond coating of the ET20D tip increases the cutting and lateral abrasion effect.



ET40 (1 per pack)

ACF88012

Long retreatment tip, 40mm, taper 4%

The ET40 is a steel tip for rapid removal of broken instruments in the middle third of wide, straight canals.



ET40D (1 per pack)

ACF88014

Long retreatment tip, 40mm, diamond-coated 30 μ m, taper 4%

The ET40D is a diamond-coated steel tip for retreatment of very hard material in the middle third.



ET25S (Short 15mm, 4% Taper) (1 per pack)

ACF88021

ET25 (Standard 20mm, 3% Taper) (1 per pack)

ET25L (Long 25mm, 3% Taper) (1 per pack)

ET25 Series Titanium-Niobium tips

ET25S - Designed for retreatment in the coronal third and the isthmuses.

ET25 - Suggested tip for retreatment in the middle and apical thirds and the extraction of broken instruments.

ET25L - The ET25L (long) is suitable for retreatment in the apical third and long, straight canals.



Apical Surgery



AS3D 3mm / 9% Taper (1 per pack)	ACF00065
AS6D 6mm / 9% Taper (1 per pack)	ACF00079
AS9D 9mm / 8% Taper (1 per pack)	ACF00067

Diamond-coated universal tips 30µm

The AS3D tip is used for apical surgery of anterior teeth. It should be used at the lowest possible effective power setting. The AS6D should be used to obtain a preparation length of at least 5mm. The AS9D tip is for preparation of the coronal third.

The AS9D tip should first be introduced into the canal and oriented in the root axis before being activated to prevent the creation of a "false route".



ASRD Right Orientated Tip (1 per pack)	ACF00081
ASLD Left Orientated Tip (1 per pack)	ACF00080

Diamond-coated 30µm, length 3mm, taper 10%

The ASRD and ASLD tips are recommended for apical surgery of premolars and molars. These instruments should be used with very light pressure.



Retro Surgery



S12-70D (1 per pack)	ACF00118
-----------------------------	----------

Retro surgery tip angled at 70°, diamond-coated 30µm, length 5mm, taper 9%

The S12-70D tip is recommended for the treatment of posterior areas, in canals that are difficult to access or roots with specific orientations.



P14D (1 per pack)	ACF00106
--------------------------	----------

Universal retro surgery tip diamond-coated 30µm, length 5mm, taper 7%

The P14D tip is recommended for the preparation of canals in anterior teeth.

The micro-retro tips make minimum treatment possible providing fast healing.



P15LD Left Oriented Tip (1 per pack)	ACF00107
P15RD Right Oriented Tip (1 per pack)	ACF00108

Retro surgery tips, diamond-coated 30µm, length 5mm, taper 7%

The P15LD and P15RD tips are recommended for the preparation of premolar and molar canals.





ENDO Success Kit

Kit Contains: ET18D, ETBD, ET20, ET25S, ET25, ETPR, Universal Wrench & Autoclavable Stand.


Endo Success Kit	ACF00737
------------------	----------



ENDO Success Apical Surgery Kit

Kit Contains: AS3D, AS6D, AS9D, ASLD, ASRD, Universal Wrench and Autoclavable Stand.

Endo Success Apical Surgery Kit	ACF00069
---------------------------------	----------



Retro Kit

Kit Contains: S12/70D, AS6D, ASLD, ASRD, 4 Yellow Wrenches and Autoclavable Box.

Retro Kit	ACF00910
-----------	----------



Micro Retro Kit

Kit Contains: P14D, P15LD, P15RD, 3 Yellow Wrenches and Autoclavable Box.

Micro Retro Kit	ACF00912
-----------------	----------



ENDO Success Canal Access Prep Kit

Kit Contains: CAP1, CAP2, CAP3, Universal Wrench and Autoclavable Stand.

Endo Success Canal Access Prep Kit	ACF88180
------------------------------------	----------



Autoclavable Stand

Autoclavable stand with silicone inserts, holds 6 tips.

Autoclavable Stand	ACN88449
--------------------	----------

Scaling & Restoration

Scaling Tips 1 per pack

#1S	Slim scaling tip	ACF00245
#1	Universal Scaling Tip	ACF00246
#2	Heavy Calculus Tip	ACF00247
#3	Stain Removal Tip	ACF00248
10P	Shallow Pockets Tip	ACF00253
10X	Interproximal Scaling Tip	ACF00359
10Z	Subgingival Scaling Tip	ACF00254

Excavus Cavity Preparation Tips 1 per pack

EX1	Ball Diamond Tip	ACF02040
EX2	Mesial Half Ball Diamond Tip	ACF02041
EX3	Distal Half Ball Diamond Tip	ACF02042
EXR	Right Half-Ball Diamond	ACF02043
EXL	Left Half-Ball Diamond Tip	ACF02044

Perfect Margin Veneers Tips 1 per pack

PMV1	Diamond Ball Tip	ACF02021
PMV2	Half Ball Diamond Ext. Spoon Tip	ACF02022
PMV3	Half Ball Diamond Int. Spoon Tip	ACF02023
PMV4	Half Ball Polishing Smooth Ext. Spoon Tip	ACF02024
PMV5	Half Ball Polishing Smooth Int. Spoon Tip	ACF02025
PMV6	Non Diamond Polishing Smooth Ball Tip	ACF02026

Perfect Margin Rounded Tips 1 per pack

PM1	1/4 Round Diamond Preparation Tip	ACF02250
PM2	1/4 Round Diamond Finishing Tip	ACF02251
PM3	1/4 Round Smooth Polishing Tip	ACF02252
PM4	1/4 Conical Diamond Corono-radicular Prep. Tip	ACF02253

Perfect Margin Shoulder Tips 1 per pack

PMS1	1/4 Round Diamond Preparation Tip	ACF02254
PMS2	1/4 Round Diamond Finishing Tip	ACF02255
PMS3	1/4 Round Smooth Polishing Tip	ACF02256

Loosening Tips 1 per pack

C20	Inlay / Onlay Tip	ACF00113
5AE	Crown Removal Tip	ACF00249
ETPR	Loosening Crowns & Posts Tip	ACF88019

Periodontal

Perio & Root Planing Tips 1 per pack

H1	Universal Diamond Perio Tip	ACF00366
H2L	Left Diamond Perio Tip	ACF00367
H2R	Right Diamond Perio Tip	ACF00368
H3	Universal Curette Tip	ACF00369
H4L	Left Curette Tip	ACF00114
H4R	Right Curette Tip	ACF00115

Biofilm Disruption Tips 1 per pack

TK1-1S	Short BDR Tip	ACF01001
TK1-1L	Long BDR Tip	ACF01004
TK2-1L	Left BDR Tip	ACF02162
TK2-1R	Right BDR Tip	ACF02161

Perio Soft Tips 4 per pack

PH1	Universal Carbon Tip	ACF00702
PH2L	Left Carbon Tip	ACF00705
PH2R	Right Carbon Tip	ACF00706

Perio Precision Tips 1 per pack

P2L	Perio Precision Left Tip	ACF00090
P2R	Perio Precision Right Tip	ACF00091

Perio Fine Tips 1 per pack

PFU	Perio Fine Tip Universal	ACF02170
PFL	Perio Fine Left	ACF02171
PFR	Perio Fine Right	ACF02172

Implant Protect Titanium Tips 1 per pack

IP1	Universal Implant Cleaning Tip	ACF02121
IP2L	Left - For use on medium implant threads	ACF02122
IP2R	Right - For use on medium implant threads	ACF02123
IP3L	Left - For use on narrow implant threads	ACF02124
IP3R	Right - For use on narrow implant threads	ACF02125

Kits

Ultrasonic Tip Kits

Scaling Kit - 1, 2, 10X, 10Z	ACF00900
Perio Kit - H3, H4L, H4R	ACF00902
Bio-Film Disruption Kit - TK1-1S, TK1-1L, TK2-1L, TK2-1R	ACF00904
Perio Soft Kit - PH1, PH2L, PH2R	ACF00906
Implant Protect Kit - IP1, IP2L, IP2R, IP3L, IP3R	ACF02120
Perfect Margin Rounded Kit - PM1, PM2, PM3, PM4	ACF00738
Perfect Margin Shoulder Kit - PMS1, PMS2, PMS3, PM4	ACF00736
Perfect Margin Veneers Kit - PMV1 to PMV6 (6 tips)	ACF02020
Excavus Kit - EX1, EX2, EX3, EXL, EXR	ACF00739



BUC Endodontic Ultrasonic Tips

Designed for Access Refinement 1 per pack



BUC 1	The Cutter	OS930-041
BUC 1A	The Slim Cutter	OS930-048
BUC 2	The Planer	OS930-042
BUC 2A	The Precision Planer	OS930-046
BUC 3	The Digger	OS930-043
BUC 3A	The Slim Digger	OS930-049

CPR Endodontic Ultrasonic Tips

Designed for Retreatment 1 per pack



CPR 1	The Mover	OS930-011
CPR 2D	The Workhorse	OS930-012
CPR 3D	The Problem Solver (Small)	OS930-013
CPR 4D	The Problem Solver (Medium)	OS930-014
CPR 5D	The Problem Solver (Large)	OS930-015
CPR 6	Violet 20mm	OS930-016
CPR 7	Blue 24mm	OS930-017
CPR 8	Bronze 27mm	OS930-018

KiS Endodontic Ultrasonic Tips

Designed for Surgery 1 per pack



KiS 1D	The General	OS931-011
KiS 1D4	The General - Longer 4mm cutting surface	OS952-614
KiS 2D	The General - For larger roots	OS931-012
KiS 3D	The Buccal (Left & Right)	OS931-013
KiS 4D	The Lingual (Left & Right)	OS931-014
KiS 5D	The Buccal (Left & Right)	OS931-015
KiS 6D	The Lingual (Left & Right)	OS931-016



Pac-Dent Endo Aid Kit

Finger held autoclavable organizer that gives dental clinicians complete control over endodontic files.

Kit Contents: 1 organizer with integrated ruler, 2 rings to fit all finger sizes, 50 disposable triangular foam refills.

Endo Aid Kit	PD802
50pk Triangular Lavender Foam Refills - Medium	PD803L
50pk Triangular Teal Foam Refills - Soft	PD803T
50pk Triangular Grey Foam Refills - Firm	PD803G
Endo Organizer (ring w/ruler)	PD803R



Jordco EndoRing II Kit

The EndoRing II organizer enables the clinician to directly place, store, measure and clean endodontic hand and rotary instruments within the operating field. It reduces the risk of cross-contamination and minimizes the passing of sharps between dentists and staff.

EndoRing II with Metal Ruler, Blue	JCER2
EndoRing II with Reduced Glare Metal Ruler, Black	JCER2GYRG
EndoRing II without Metal Ruler, 3 color, Variety Pack	JCER2CVP
Jordco e-Foam, Blue/Yellow Combo, 48pk	JCERF
Jordco e-Foam, Blue, 48pk	JCERFB
Jordco e-Foam, Yellow, 48pk	JCERFY
Jordco e-Foam, Silver, 48pk - HIGH DENSITY	JCERFHD



Jordco FileCaddy

The FileCaddy provides a unique bulk storage solution for all endodontic files. The clinician places multiple files securely into the special foam insert, closes the lid and steam autoclaves the entire organizer.

Each kit contains: 1 EndoRing FileCaddy, 12 Foam Insert Refills



Jordco FileCaddy	JCEFC
Jordco FileCaddy Foam Inserts 12pk	JCEFCF



Endo Sponges

Fully disposable blue polyurethane sponges. Each package contains 50 pieces of sponge.

Round Endo Sponge (2" dia. X 3/8" thick)	PD811
Round Endo Sponge (2 1/2" dia. X 3/8" thick)	PD812
Square Endo Sponge (2 1/8" X 2 1/8" X 1/4" thick)	PD813



iMatrix Sectional Matrix System

The iMatrix Sectional Matrix System aids in the production of anatomically correct Class II composite restorations. It features ultra-retentive NiTi matrix rings that deliver non-slipping contact and exert ideal tooth separation force.

The iMatrix™ Sectional Matrix System offers a starter kit (with instruments), intro kit (without instruments for those that already have instruments) and refill kit configurations.

iMatrix Starter Kit Contents: x1 Standard Ring, x1 Narrow Ring, x35 each of Small, Medium & Large Wedges, x35 each of 3.5, 4.5, 5.5 & 6.5mm matrices, x1 Forceps, x1 Tweezers.

iMatrix Starter Kit	PDPM01
iMatrix Intro Kit (No Instruments)	PDPM02
Interproximal Wedges - Small 100pk	PDPMWS
Interproximal Wedges - Medium 100pk	PDPMMW
Interproximal Wedges - Large 100pk	PDPMWL
3.5mm Matrices Refill 100pk	PDPMM3
4.5mm Matrices Refill 100pk	PDPMM4
5.5mm Matrices Refill 100pk	PDPMM5
6.5mm Matrices Refill 100pk	PDPMM6
Ring Refill 2pk - 1 Standard (Black) & 1 Narrow (Grey) Ring	PDPMR1
Ring Refill 2pk - 2 Standard (Black) Rings	PDPMR2
Ring Refill 2pk - 2 Narrow (Grey) Rings	PDPMR3
iMatrix Tweezers	PDPMT1
iMatrix Ring Forceps	PDPMF1

Pulpdent Rainbow Wooden Wedges

Strong wooden wedges in four colour coded sizes. Anatomically contoured to avoid laceration of gingival tissue. Available in an assorted starter kit or in individual refill packs.

Starter Kit (x100 of each size/color), 400pk	PUW
Blue Wedge Refill Pack (12mm) 100pk	PUWB
Green Wedge Refill Pack (13mm) 100pk	PUWG
Yellow Wedge Refill Pack (15mm) 100pk	PUWY
White Wedge Refill Pack (17mm) 100pk	PUWW



Plug In Bur Stands

The silicone plugs hold any type of bur securely. Made from stainless steel and silicone. Can be used with all known decontamination, cleaning and sterilization processes. Each stand can be rotated to offer the burs flat or at an angle. Will hold Friction Grip or Right Angle Slow Speed burs no matter their diameter. Made in France.

Plug In Bur Stand - 6 Holes - Blue or Yellow Inserts	NM400341B/Y
Plug In Bur Stand - 8 Holes - Blue or Yellow Inserts	NM400342B/Y
Plug In Bur Stand - 12 Holes - Blue or Yellow Inserts	NM400343B/Y
Plug In Bur Stand - 16 Holes - Blue or Yellow Inserts	NM400344B/Y
Plug In Bur Stand - 20 Holes - Multi Color Inserts	NM400345
Plug In Bur Stand - 30 Holes - Multi Color Inserts	NM400346

Plug In Bur Stand - 6 Holes - 40mm Height - Blue or Yellow Inserts	NM400348B/Y
Plug In Bur Stand - 6 Holes - 50mm Height - Blue or Yellow Inserts	NM400347B/Y

Plug In Blue Silicone Inserts - 1ea	NM190200B
Plug In Yellow Silicone Inserts - 1ea	NM190200Y
Plug In Red Silicone Inserts - 1ea	NM190200R
Plug In Green Silicone Inserts - 1ea	NM190200G



12 Hole Endo Ultralight Cassette

A very strong and light endo holder, when closed the Ultralight endo holder becomes a cassette giving a high level of protection for your endodontic instruments. Electro polished to give a bright shiny surface finish. Made of stainless steel and silicone. Made in France.

12 Hole Endodontic Ultralight Cassette	NM400376
--	----------



32 Hole Endo Cassette

A stainless steel perforated box with a teflon insert, 1 endometer and a self locking cover system. Ideal for use in ultrasonic baths, cleaners and disinfectors. The box includes a silicone cursor on the cover for counting uses. Size 6.5 x 4.3 x 5cm. Made in France.

32 Hole Endo Cassette	NM182150
-----------------------	----------



Friction Grip & Right Angle Bur Blocks

Fully autoclavable aluminum blocks. Available in 18 (All FG) and 27 (18FG + 9RA) bur block styles.

18 Hole Bur Block - (18FG)	PD071
27 Hole Bur Block (18FG + 9RA)	PD072



Gelato Topical Anesthetics

Gelato Topical Anesthetic Gel is a 1oz 20% benzocaine anesthetic gel that is fast acting with no systemic absorption.

Cherry	KYGC
Grape	KYGG
Mint	KYGM



Vital-Ice Pulp Vitality Refrigerant Spray

Vital-Ice has a mild peppermint scent and a minimum temperature of -26.2°C. Each can is 6oz in size. Made in the USA.

Vital-Ice	PDVI
-----------	------



MicroPro Micro Applicators

Disposable, all-purpose micro applicators made with non-linting fiber ball ends. Each package contains 400 applicators.

Standard, Blue	PD901
Fine, Green	PD902
Ultra Fine, Purple	PD903
Tapered, White	PD904

MicroBrush Micro Applicators

With a stiff head designed for scrubbing and a strong bending portion for accurate placement of solution in hard to reach areas, Microbrush applicators offer dentists unparalleled convenience, consistency, and control. The non-absorbent fiber tips hold solution in suspension until applied to ensure maximum control and precision in every restoration.

Regular, 400Pk, Asst Colours	MBMRA400
Fine, 400Pk, Asst Colours	MBMFA400
Super Fine, 400Pk, Asst Colours	MBMSF400
Ultra Fine, 400Pk, Orange	MBMUO400


Cotton Pellets

Richmond Dental 100% Cotton Pellets are designed to provide outstanding quality and absorbency in a variety of uses.



Product	Description	Packaging	Code
Size 1	9/32" Diameter	2000 pellets per box	RD1
Size 2	7/32" Diameter	2550 pellets per box	RD2
Size 3	5/32" Diameter	2500 pellets per box	RD3
Size 4	1/8" Diameter	3000 pellets per box	RD4





ULTRA Bright Micro Mirrors

Available in 3 different sizes: 3mm, 5mm and 3x6mm these ULTRA bright Simple Stem micro mirror heads are ideal for endodontic surgery and implant work. Made in Germany.

ULTRA 3mm Micro Mirror Head (1Pk)	HK6181
ULTRA 5mm Micro Mirror Head (1Pk)	HK6182
ULTRA 3x6mm Micro Mirror Head (1Pk)	HK6183
Satin, Hexagonal Hollow, Lightweight 16g Mirror Handle	HK476



ULTRA Bright Black MICROflex Micro Mirrors

MICROflex micro mirrors feature a bendable stem for precise alignment. The DLC (Diamond Like Carbon) coating is antireflective to reduce visual fatigue. Made in Germany.

ULTRA MICROflex Oval, Black, Simple Stem	HKB5183
ULTRA MICROflex Teardrop, Black, Simple Stem	HKB5184
Micro Mirror Handle, ERGOtouch, Simple Stem, Black	HKB442



ULTRA Bright Front Surface Mirror Heads

The brightest and sharpest mirrors in the world. Featuring the highest reflection rate with crystal clear brightness and superior colour accuracy, without any discoloration. Made in Germany. **6 per pack (Cone Socket)**

ULTRA #4 FS Mirror	HK6165
ULTRA #5 FS Mirror	HK6166
ULTRA #4 FS Bifurcated Stem Mirror	HK6695
ULTRA #5 FS Bifurcated Stem Mirror	HK6696



ULTRA Bright Single Piece Dental Mirrors

Front Surface single piece mirrors are light weight, ergonomic and feature a matte smooth finish. Each instrument is autoclavable, acid resistant and suitable for thermal disinfectors. Available in #4 and #5 sizes. Made in Germany.


ULTRA #4 Single Piece Dental Mirror, Grey (10Pk)	HK7306
ULTRA #5 Single Piece Dental Mirror, Grey (10Pk)	HK7316
ULTRA #4 Single Piece Dental Mirror, Pink (10Pk)	HK7307
ULTRA #5 Single Piece Dental Mirror, Pink (10Pk)	HK7317

High Quality Mirror Handles

Each mirror handle is fine polished and all edges are smoothly rounded for comfort. Made in Germany.

1 per pack (Cone Socket)

Stainless Steel Polished, Octagonal, 24g	HK415
Stainless Steel Polished, Octagonal, Serrated, 24g	HK421
Stainless Steel Polished, Round, Serrated, 26g	HK427



Surgical Aspirator Tips

For removing zinc oxide dental cements, impression materials and waxes from vinyl furniture, face and hands, and whenever a strong, organic cleaner is required.

Large Aspirator Tips, Green, 25pk	PD008
Standard Aspirator Tips, White, 25pk	PD009
Small Aspirator Tips, Blue, 25pk	PD010



Armor Air/Water Syringe Tip & Sleeve

Armor minimizes preparation time by combining the air/water syringe tip and sleeve placements into a one-step operation. It features the acclaimed disposable syringe tip design with 6 air-core design, and a non-slip sleeve made with special copolymer film that provides superior infection control.

Armor Air Water Syringe Tip & Sleeve, 144pk	PD101-144
---	-----------



Mini Brush Applicators

Disposable mini brush applicators with a bendable tip available in 4 colors. **100 per pack**

Blue Applicators	PD500MB
Green Applicators	PD500MG
Red Applicators	PD500MR
Yellow Applicators	PD500MY



Digital X-Ray Sensor Sheaths

These specially designed sheaths with peel off paper backing are custom fitted for each digital X-Ray sensor type. The sheaths' films are made of a resin blend that is very durable. 100% latex free. **500 per pack**

MODEL	SIZE	COMPARES TO TIDI #	CODE
Gendex, XDR	Size 1	-	PDX405
Gendex, XDR	Size 2	-	PDX406
Dexis	-	20999	PDX999
Kodak 6100	Size 0	20977	PDX977
Kodak 6100	Size 1	20978	PDX978
Kodak 6100	Size 2	20979	PDX979
Sirona	Size 2	20890	PDX890
Sirona	Size 1	20904	PDX904
Schick	Size 1	20824	PDX824
Schick	Size 2	20825	PDX825
Suni/Lightyear	Size 2	20819	PDX819



NeoTip Air/Water Syringe Tips

Comparable to Crystal Tips®, the inner tube does not move under pressure. A groove in the tip ensures a snug locking fit so that the tip stays in place throughout the procedure.

NeoTip Air/Water Syr Tip Asst 250/pk	PD105-250AST
NeoTip Air/Water Syr Tip Asst 1500/pk	PD105-1500AST



Accessories

Gutta Cutter	GC
Tooth Slooth, 4 per pack	TS
Tooth Slooth Blue II, 4 per pack	TSB
Wonder Orange Solvent, 8oz Bottle	PUWO
DiaDent Assorted Silicone Rubber Stops, 120 per pack	D701-291
DiaDent Dia-Gauge	D781-001

Endo/Tech proudly seeks opportunities to help our world and our professional community. We contribute product to organizations who give back in their own communities, both near and far. We recycle and compost and have reduced our paper usage through electronic invoicing. Endo/Tech makes annual donations on behalf of each staff person, to a charity of their choosing, and we support the future stars of dentistry by donating and contributing to dental schools and dental education.



389 - 5991 Spring Garden Road
Halifax, Nova Scotia B3H 1Y6
Canada

Toll-Free Phone: 888.554.3636

Toll-Free Fax: 877.553.8329

www.endo-tech.com / endotech@endo-tech.com



Find us on Twitter
[@EndoTechRGS](https://twitter.com/EndoTechRGS)



Endo/Tech reserves the right to modify products and prices without notice. Products may not be shown to scale.
Printed in Canada. Copyright 2023, Endo/Tech. All rights reserved.

The product names used in this catalog are for identification purposes only. All trademarks and registered trademarks are the property of their respective owners.